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Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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Important information for users of KARL STORZ devices

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1 Important information for users of KARL STORZ equipment

Thank you for your confidence in the brand name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the HYSTEROMAT® E.A.S.I.. All required details and all actions on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

General warning

Users of the HYSTEROMAT® E.A.S.I. must be thoroughly trained in its use for the applicable medical procedures. The device and accessories should not be used until you have familiarized yourself fully with the instruction manual, assembly, and the care instructions. Instruction manuals should be carefully studied and be available to the surgical team during the procedure. It is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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

E.A.S.I. =
Endoscopic Automatic System for Irrigation

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2 Example of a workstation
(gynecology)

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




*SCB = STORZ Connection Bus

**4 Controls, displays,
connectors and their uses**

- ❶ Power switch, illuminated
0 = off
- ❷ Display and touch screen
- ❸ Irrigation pump with locking lever
- ❹ Suction pump with locking lever
- ❺ Pump roller
- ❻ Pressure sensor for the pressure measuring chamber of the irrigation tube
- ❼ Tube holder for the suction tube
- ❽ Footswitch connection socket
- ❾ Connection for shaver
- ❿ Connections for SCB*
- ⓫ Potential equalization connector
- ⓬ Power cord socket
- ⓭ Fuse holder

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Read the instructions carefully before operating the
equipment!

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Potential equalization connector

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Applied part type BF in accordance with
IEC 60601-1

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Electronic information product pollution control
(China RoHS)

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This device has been labelled in accordance with the
European Directive 2002/96/EC on waste
electrical and electronic equipment (WEEE).

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This material is intended for single use and must
not be resterilized.

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1	Important information for users of KARL STORZ equipment	4
2	Example of a workstation (gynecology)	5
3	Photographs of the equipment	6
4	Controls, displays, connectors and their uses	7
5	Symbols employed	9
6	General information	13
6.1	The HYSTEROMAT® E.A.S.I. device	13
6.2	Intended use	13
6.3	User qualification	14
6.4	Operating modes	14
6.5	Operation	14
6.6	Key performance features	14
7	Warnings and cautions	15
7.1	Definitions	15
7.2	Training in the operation and function of the device	21
7.3	Patient profile	22
7.4	User profile (specialist physician)	22
7.5	User profile (assistant)	23
7.6	Use	23
7.7	User position	24
7.8	Safety precautions at the site of installation	25
7.9	Safety precautions when operating the device	25
8	Explanation of terms (in alphabetical order)	26
9	Installation and operating instructions	31
9.1	Unpacking the equipment	31
9.2	Basic equipment	31
9.3	Installing and connecting up the device (preparations)	31
9.4	Operating the unit	34
9.5	Preparing connection of the irrigation tube	35
9.6	Preparing connection of the suction tube	38
10	Menu description	41
10.1	Components of the header/footer	41
10.1.1	Components of the display	41
10.2	Start screen	42

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6 General information

6.1 The HYSTEROMAT® E.A.S.I. device

The HYSTEROMAT® E.A.S.I. is a combined pressurized suction and irrigation pump with a continuous fluid flow for hysteroscopic interventions. Alternatively, the two pumps can be operated independently of each other with flow regulation, e.g. for laparoscopic interventions.

In pressurized mode, the control compensates for the flow-related fall in pressure between the supply and the opening of the hysteroscope both vis-a-vis the display and setpoint. This allows the user to work directly with intra-uterine pressure.

The device accumulates the irrigation liquid supplied and can inform the user when preset limit values have been exceeded (monitoring of irrigation fluid supply).

The HYSTEROMAT® E.A.S.I. is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface. The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

6.2 Intended use

Suction/irrigation pumps and their accessories are used for the introduction of irrigation fluids into organs and operating fields as well as the suctioning off of irrigation and bodily fluids, secretions and tissue during endoscopic interventions.

For safety reasons, the device must only be used as intended. Use of the HYSTEROMAT® E.A.S.I. is contraindicated if, in the opinion of the responsible physician, the patient's health is endangered through its use.

The HYSTEROMAT® E.A.S.I. may only be operated with tubing sets, sheaths and accessories which have been designated as suitable for the device by KARL STORZ.

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General information

This applies in particular for tubes and sheaths, since the device uses the predetermined characteristics of these elements for determining the intra-uterine pressure. No responsibility for the safe operation of the device can be accepted if any accessories other than those specified are used.

The device can be used in combination with a HF surgical device.

6.3 User qualification

The HYSTEROMAT® E.A.S.I. may only be used by physicians and trained medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the device.

6.4 Operating modes

The HYSTEROMAT® E.A.S.I. has 3 operating modes for surgical and diagnostic interventions/operations. The operating modes are 'HYS' (hysteroscopy) for gynecological diagnostics, 'RES' ('resection') for surgical hysteroscopic interventions and 'LAP' ('laparoscopy') for laparoscopic interventions.

6.5 Operation

The device features a color monitor with a touch-sensitive screen for making entries (touch screen). The information displayed and the possible entries are adapted to the specific work flow status. Both pump beds feature a rotating tensioning lever with which the user can open or close the pump bed for insertion and removal of the tubes

6.6 Key performance features

In the case of first faults, pressure specification, display and regulation function as per the technical specification in the instruction manual OR an alarm is issued.

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7 Warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this instruction manual, they should be carefully reviewed to ensure the safe and effective operation of this device. To make the signal words stand out more clearly, they are accompanied by a pictogram.

7.1 Definitions

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: Notes contain additional information on the safe operation of the device.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel or yourself at risk.

WARNING: The electrical installations in the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep the unit out of reach of patients.

WARNING: For safety reasons, when the system is in use the output sockets of the device and the patient must not be touched at the same time.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

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Safety instructions

Warnings and cautions

WARNING: Combinations of medical devices are only assured to be safe if:

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

WARNING: To avoid the risk of an electric shock, this unit may only be connected to a power supply network with a protective conductor.

WARNING: When used in combination with HF surgical devices: Only HF operating modes with a maximum output voltage of 5.5 kVpp may be selected.

WARNING: Grounding reliability can only be achieved when the equipment is connected to 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: Test this device prior to each surgical procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

WARNING: Always unplug the device before performing any maintenance work on it.

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WARNING: Do not open the device! Danger of an electrical shock. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only (see §4 German Medical Device Directive). Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims.

WARNING: Use of electromechanical medical devices can cause electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions must be taken to ensure that the patient's wellbeing is maintained in the event of such interference.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

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Safety instructions

Warnings and cautions

WARNING: Reel pump (pressure-controlled pump motor). Do not touch the pump reels 5 when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

WARNING: The irrigation tubing set must be filled with irrigation fluid before use on a patient.

WARNING: The fluid level in the irrigation bag must be monitored carefully to ensure that no air is pumped into the cavum. Another bag of irrigant should be started prior to emptying the bag being used.

WARNING: The intracavitary pressure is always equal to or lower than the pressure at the outlet of the HYSTEROMAT® E.A.S.I. which is used to regulate the device.

WARNING: The user should always use the lowest possible pressure for sufficient distention and/or irrigation.

WARNING: Do not reuse disposable tubing sets. Discard after use.

WARNING: If a footswitch no longer reacts when pressed, the device must not be operated.

WARNING: Before you activate the pump, make sure that the tubing sets have been inserted properly and that the pressure measuring chamber is correctly positioned.

WARNING: The time display showing the seconds in the header indicates that the device control is working properly. If the changing seconds are not visible, the device must not be operated.

WARNING: Before you start the intervention, make sure that the tubing sets used are not leaking or kinked and that the ends of the tubes are properly secured.

WARNING: The device regulates based on the pressure measured at the outlet of the irrigation pump. The pressure in the cavity depends, among other things, on the height difference between the device and cavum. Try to keep the difference as small as possible.

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Safety instructions

Warnings and cautions

WARNING: Any treatment may only be performed if visual observation of the action of the device is ensured.

WARNING: Attaching a pressurized measuring chamber (cambered diaphragm) can result in an erroneous pressure measurement. The irrigation tube should always be unpressurized when installed.

WARNING: Observe the applicable local regulations when disposing of disposable tubing.

WARNING: Observe the availability signal during power-up. If it can't be heard the device must not be used.

WARNING: Screw the patient tube's connector firmly into the pump tube's coupling to avoid leakages.

WARNING: Servicing may only be performed by specialists authorized by KARL STORZ.

WARNING: The device displays the intra-uterine pressure. The correct functionality depends on correct selection of the used sheath. Make sure to use matching sheath and procedure.

WARNING: The HYSTEROMAT® E.A.S.I. may be subject to sudden failures. Therefore you should ensure that a replacement unit is available or be prepared to apply alternative procedures if the intervention requires permanent fluid transportation/distension.

WARNING: The HYSTEROMAT® E.A.S.I. relies on the characteristics of the used sheath for correct pressure measurement. Never use the device with sheaths other than specified in the instruction manual.

WARNING: The one-day tubing set for irrigation can remain in the reel pump for a maximum of one surgery day (24 hours) during which it can be used for several operations. Afterwards it must be disposed of.

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Safety instructions

Warnings and cautions

WARNING: Valves and tube clamps should be either completely opened or completely closed.

WARNING: Make sure that the set volume of the loudspeakers allows for perception of acoustic signals over normal background noise.

WARNING: If the touch screen no longer reacts when pressed, the device must not be operated.

WARNING: The device must not be operated if the touch screen or display are defective.

WARNING: It is not permitted to press more than one touch point simultaneously.

WARNING: The device may only be operated under the environmental conditions specified.

WARNING: Please ensure that there is sufficient air circulation for your device.

WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: The device may only be operated at the voltage stated on the rating plate.

WARNING: When cleaning the device, always observe the instructions in the instruction manual (see section "Cleaning and care").

WARNING: When replacing fuses, only use fuses of the correct rating.

WARNING: The Bipolar Instrument Mode is designed for the compensation of the RF instruments being used. Do not use while no RF instruments are in use.

CAUTION: Avoid allowing liquids to enter the housing. Do not store liquids on or above the unit. If, despite precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies to the formation of condensation).

CAUTION: When cleaning the device surface, the area around the pressure sensor ⑥ should be wiped carefully due to its sensitive surface.

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Safety instructions

Warnings and cautions

NOTE: Observe the applicable local regulations when disposing of disposable tubing.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

NOTE: This unit may only be operated with the Irrigation Tubing Set 031717-01 intended for this purpose and which is indicated by the bar on the packaging.

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Before using the unit on the patient, it is imperative that you be acquainted with how the unit operates and is controlled.

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7.2 Training in the operation and function of the device

Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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7.3 Patient profile

Species: Human
Ethnic group, race, skin color: irrelevant
Nationality, language: irrelevant
Gender:
LAP: irrelevant; HYS: female; RES: female
Age: irrelevant
Weight: age appropriate
Health: suitable for endoscopic interventions
Mental health:
irrelevant, patient is never the user
Affected body part or tissue type
HYS:
Body part: uterus
Access port: vagina
State: suitable for hysteroscopic interventions
RES:
Body part: uterus
Access port: vagina
State: suitable for hysteroscopic interventions
LAP:
Body part: abdominal cavity
Access port: abdominal wall
State: suitable for laparoscopic interventions

7.4 User profile (specialist physician)

- Recognized medical skills in the area of gynecology/ laparoscopy (specialist physician).
- Adequate powers of comprehension to assess the situation rationally at all times during the application.
- Adequate language skills in at least one of the languages used on the device and in the instruction manual.
- Be thoroughly trained in the operation and use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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7.5 User profile (assistant)

- Completed training as medical assistant.
- Adequate language skills in at least one of the languages used on the device and in the instruction manual.
- Be thoroughly trained in the operation and use of the device.
- Have been trained thoroughly at least once in the use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

7.6 Use

The device may only be used in hospitals and medical practices in the ambient conditions given in the technical data.

- Not on intensive care units.
- Not in mobile facilities (cars or trucks, ambulances, aircraft, ships or trains).
- Position it on a vibration-free, fixed base.
- Visual distance from the physician to the device when operated: max. 2 m.
- Visual distance for operating the device: 30-70 cm.
- Viewing angle to the front of the device: max. $\pm 45^\circ$.
- The use of the device requires the division of work.
- Frequency of use: up to several times a day.
- Length of use: from a few minutes to several hours a day.
- Can only be moved by installing it on an equipment trolley.
- Can only be swiveled by installing it on ceiling supply units.

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Safety instructions

- Patient must not necessarily be under general anesthesia during the intervention (local anesthesia).
- May be used with other devices (light source, HF surgery etc.).
- To be operated in ambient light between 100 Lux and 1500 Lux.
- Ambient temperature in immediate proximity of the device during operation: max. 40°C.
- Sound pressure from background noise < 60 dBA.
- The device is placed outside the sterile area.
- Must only be operated with parts for single use.

7.7 User position

- For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$ at a distance of approx. 30 – 70 cm to the front panel.
- For observation of the actual values and for operation of the footswitch during the application, a visual distance to the device of 2 m is assumed (the tube length is 2.5 m).

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
7.8 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The device is equipped with a connector for potential equalization. It should be connected up in accordance with the applicable national regulations.

WARNING: The device must not be used in oxygenated environments. It is not intended for use in hazardous zones.

This means, among other things:

WARNING: When using easily combustible and explosive inhalation anesthetics or mixtures thereof, the device may not be operated in the demarcated hazard zone.

This also applies to easily combustible and explosive chemicals, e.g. skin disinfectants and surface disinfectants.

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7.9 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operates properly before use. During treatment with the HYSTEROMAT® E.A.S.I., the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring vital signs and the anesthetic.

Any treatment may only be performed if visual observation of the action of the device is ensured.

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8 Explanation of terms (in alphabetical order)

И **NOTE:** The terms are in alphabetical order in all languages. As such, it is possible that the German text on one page does not correspond to the English or Spanish text in terms of content.

Activated: When the line voltage is switched on, the user can move the rollers of the irrigation and suction pump via a button on the screen. This is referred to as 'activation of the pumps'.

И **NOTE:** Further preconditions must be satisfied prior to activation. These are described as of page 31 (Installation and operating instructions).

Bipolar Instrument Mode: The Bipolar Instrument Mode is selected by the user if an increased irrigation pressure is necessary, e.g. if an HF electrode is inserted into the working channel of the instrument (endoscope), resulting in a reduction of the instrument's cross section for irrigation liquid.

! **WARNING:** The Bipolar Instrument Mode is designed for the compensation of the RF instruments being used. Do not use while no RF instruments are in use.

Blood/Mucus Mode: In Blood/Mucus Mode, the flow is increased to enable improved suction of blood and mucous particles, thus providing the user with a better view. The Blood/Mucus Mode can be set in the operating menu 'Pressure/Flow Settings' (only with Continuous Flow procedures). In the Blood/Mucus Mode, the pressure can be increased in percentages by the user. An increase in pressure is required should it be necessary to prevent the cavum from collapsing or in order to push back or stop bleeding. If, for example, the user has 60 mmHg and increases this by 50%, this is then 90 mmHg.

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Continuous Flow: Continuous Flow CF: In the Continuous Flow Mode the pump firstly irrigates in order to build up pressure. The suction pump is inactive until the pressure has been achieved. Once the pressure in the uterus has been achieved, the Continuous Flow Mode is activated. Both pumps then irrigate and suck with the same flow rate (e.g. 50 ml/min). As such, there is no fall in pressure.

The Continuous Flow Mode offers the advantage that clear fluid is continuously conveyed. In contrast to this, with the Single Flow Mode the fluid in the uterus tends to be cloudy or bloody.

Day set / day tube set: The day set concerns only the irrigation tube and not the suction tube.

Purpose of the day set: There is an interface in the irrigation tube set; this is kept sterile. Before each intervention a new patient tube (length approx. 2 m) is connected to this interface. After each operation, rather than replacing the entire irrigation tube only the patient tube of the irrigation tube set is replaced.

This means that the costly sterile parts of the irrigation tube set (e.g. pressure measuring chamber) only have to be replaced once a day; only the non-sterile (e.g. contaminated with blood) patient tube is replaced after every intervention.

Deactivated: Opposite of → activated.

Double Flow Mode: This mode is available for laparoscopy. Here, the suction and irrigation flow rates can be changed and the respective pumps activated and deactivated all independently of each other.

Flow modes: There are three different flow modes:

- → Continuous Flow Mode
- → Single Flow Mode
- → Double Flow Mode (for laparoscopy)

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Explanation of terms

in alphabetical order

LAP Mode: The LAP procedures are exclusively suction and irrigation procedures.

A pneumoperitoneum for laparoscopic procedures is created via insufflation, e.g., with carbon dioxide gas.

Menus: A distinction is made between

- Operating menus (the user works with these during the intervention – e.g. Normal Mode operating menu)
- 'Settings' menu: this is used to set all basic settings (e.g. prior to the intervention), e.g. pressure/flow settings, monitor language settings etc.

Normal Operation: During Normal Operation, regulation is based on the intra-uterine pressure. This is done using characteristics which are saved in the device for the specific sheaths. These characteristics differ depending on the type of intervention (hysteroscopy or resection) and the associated sheath.

The intra-uterine pressure is calculated automatically in Normal Operation. This means that the user does not have to worry about falls in pressure from the irrigation liquid bag via the pressure measuring chamber to the instrument. The pressure shown on the monitor in the operating menu 'Normal Operation' is the actual intra-uterine pressure.

In the operating menu 'Normal Operation', the user can only change the pressure and not the flow. However, using the menu 'Pressure/Flow Settings', he can set a specific flow for a specific pressure. The flow is then automatically supplied by the system during irrigation depending on the pressure selected.

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Explanation of terms

in alphabetical order

Operating state: There are three different pump operating states for hysteroscopy and resection interventions.

- Normal Operation (when neither the Bipolar Instrument Mode nor the Blood/Mucus Modus has been activated)

- → 'Bipolar Instrument' Mode

- → Blood/Mucus Mode

Only the Normal Operation Mode is available for laparoscopic interventions.

Pressure increase: In the → Blood/Mucus Mode, the pressure can be increased in percentages by the user. An increase in pressure is required should it be necessary to prevent the cavum from collapsing or in order to push back or stop bleeding. If, for example, the user has 60 mmHg and increases this by 50%, this is then 90 mmHg.

Procedures: With the Hysteromat E.A.S.I., the user works exclusively with procedures. On delivery, the device features predefined procedures for the Gyn Mode (HYS and RES) and the LAP Mode. The user cannot change, rename or delete these predefined procedures.

He can, however, copy a predefined procedure and save this under another name and then change the parameters individually (e.g. the desired pressure or flow).

The following applies: a GYN Mode procedure (HYS or RES) with Continuous Flow sheath can only be changed into a GYN Mode procedure with the same Continuous Flow sheath.

A GYN Mode procedure (HYS) with Single Flow sheath can only be changed into a GYN Mode procedure with the same Single Flow sheath.

A LAP Mode procedure can only be changed into an LAP Mode procedure.

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9.1 Unpacking the equipment

Carefully unpack the HYSTEROMAT® E.A.S.I. and its accessories from their packing. Check the delivery for missing items and possible damage. Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

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9.2 Basic equipment

- 1 HYSTEROMAT® E.A.S.I. 26340020-1
- 1 Power cord 400 A (400 B)
- 1 SCB connecting cord 20090170
- 1 Instruction manual 96126014D
- 1 Basic Set Tubing Sets 031917-10 consists of 5x Irrigation Tubing Sets 031717-01 and 5 Suction Tubing Sets 031217-01

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9.3 Installing and connecting up the device (preparations)

WARNING: Observe the section 'Safety precautions at the site of installation' on page 25.

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

Potential equalization

Connect the potential equalization cable to the potential equalization connector 9.


WARNING: The device's ground line should be installed by a qualified electrician.

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Connecting the power cord

Connect the power cord when the device is switched off. Insert the power cord into the socket 10 as far as it will go.

The device may only be operated at the voltage stated on the rating plate.

WARNING: The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

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WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.


NOTE: Ensure that access to the power plug is ensured at the site of installation.

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Connecting the footswitch

Insert the connecting cable of the footswitch into the corresponding socket 7.


NOTE: The footswitch is an optional accessory.

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Connecting the shaver

Insert the connecting cable of the shaver into the corresponding socket 9.

NOTE: The shaver is an optional accessory.

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
Installation and operating instructions

Connecting up the SCB cable (optional)

In order to use this function, insert one end of the SCB cable (both ends are identical) into the connection socket ⑧ of the HYSTEROMAT® E.A.S.I.. Connect the other end to an available SCB connection socket on a device in the existing SCB network or a device with which you wish to establish a network.

NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑧. Connect the other end of the cable to the KARL STORZ SCB control device (KARL STORZ Communication Bus) or other SCB devices (see 'SCB System' instruction manual).

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
9.4 Operating the unit

CAUTION: Before switching the HYSTEROMAT® E.A.S.I. on, make sure that no tubing set has been inserted into the device and that the cables are properly connected.

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Press the power switch ①.

! **CAUTION:** Do not touch the touch screen while you switch the device on.

Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo. After this, a white screen appears.

! **WARNING:** Check the white screen for color errors. In the event of any errors the device must not be used.

During this period, the device runs a self test. Once the test has been successfully completed, an acoustic availability signal sounds. An animation appears on the screen which assists you in inserting the tubing set.

If the self test fails, the message 'Safe State' is displayed and an alarm sounds. Verify that the device is properly connected and that the preparations have been made according to the instructions.

For further information about troubleshooting please refer to the section 'Technical description: Information and alarm signals' from page 85.

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Animation

The animation shows the order in which the tubing set should be inserted. The animation runs in an endless loop. This menu is exited when the device detects that:

- the irrigation and suction pump levers have been opened and closed and the pressure measuring chamber on the sensor ⑥ has been adapted (*Continuous Flow Mode*).
- the irrigation pump lever has been opened and closed and the pressure measuring chamber on the sensor ⑥ has been adapted (*Single Flow Mode*).

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9.5 Preparing connection of the irrigation tube

! **WARNING:** The packaging for the tubing set must not be damaged (sterility!). Check the tubing set for damage etc. before use.

Remove the irrigation tube (blue, with pressure measuring chamber) from the sterile packaging and connect it.

Pass the irrigation tube on for installation in the pump while observing aseptic procedures (safety precautions).

! **WARNING:** The connection of the irrigation tubing set (sheath) to the patient must remain in the sterile area.

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Unlocking the irrigation pump

WARNING: Do not touch the pump reels **5** when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.


Unlock the irrigation pump by turning the lever clockwise as far as it will go. As soon as the lever is actuated, the pump can no longer be activated.

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
Positioning the pressure measuring chamber

Position the pressure measuring chamber and turn it clockwise until fixed.

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Inserting the irrigation tube

Insert the irrigation tube into the pump bed.


WARNING: Make sure that the tube is completely inserted into the pump bed.

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Locking the irrigation pump

Lock the irrigation pump by turning the lever counterclockwise until it clicks into place.

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
Connect the irrigation tube to the irrigation bag(s). Make sure that no liquids are stored on or above the device (see the section 'Safety instructions', p. 15)


Note: The tubes for the irrigation bag are fitted with a blue clamp which can be used to pinch off the tubes if necessary.

! WARNING: The irrigation pump is not self-priming. As such, the irrigation bag must always be located to the side and above the device.

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Day tube sets

If the device is approved for use with day tube sets, the steps for inserting the irrigation tube set 031262-01 must be performed as described on page 35-36. The separately packed patient tube 031261-01 is first connected with the sheath and then with the irrigation tube via the coupler.

! WARNUNG: The patient tube must be disposed of after each intervention, while the tube inserted in the device may be reused.

! WARNUNG: A new patient tube must be used for each patient.

! WARNUNG: To guarantee the sterility of the connected tube set during patient changeover, the patient tubes are equipped with a separately packed, sterile plug. This must be screwed into the tube coupler from removal of the patient tube from the tube set until the next patient tube is connected.

Handling of the day tube set:

- Following the end of the intervention, the patient tube must be removed from the irrigation tube and disposed of correctly.

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- The plug is then removed from the small sterile packaging and used to seal off the open end of the pump tube.
- The new patient tube is removed from the tube set for the next patient.
- The plug in the pump's irrigation tube is now removed and disposed of, and the new patient tube is connected to the open pump tube under aseptic conditions (safety measures).
- The new plug in its sterile packaging is then offered from the sterile area.

9.6 Preparing connection of the suction tube

WARNING: The packaging for the tubing set must not be damaged (sterility!). Check the tubing set for damage etc. before use.

Remove the suction tube (red) from the sterile packaging and connect this to the instrument's sheath.

Pass the suction tube on for installation in the pump while observing aseptic procedures (safety precautions).

Warning: The connection of the suction tubing set to the patient must remain in the sterile area.

NOTE: When using the pump in Single Flow Mode, the steps for connection of the suction tube are not required. If relevant, observe the instruction manual of the suction system used.

NOTE: The tube to the sheath is fitted with a red clamp which can be used to pinch it off if necessary.

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Unlocking the irrigation pump


! **WARNING:** Do not touch the pump reels **5** when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

Unlock the suction pump by turning the lever clockwise as far as it will go. As soon as the lever is actuated, the pump can no longer be activated.

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Inserting the suction tube


Insert the suction tube into the pump bed. Make sure that the tube is completely inserted into the pump bed.

The suction tube is equipped with two tubes at the patient end which are opened using the tube squeeze valve for the respective operating mode. The tube connector at the patient end should be inserted into the suction tube bed at the top in such a way that the locking lever can be turned counterclockwise to the end position without meeting obstacles. The outflow end of the suction tube is inserted in the tube bed at the bottom.

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Locking the suction pump

Lock the suction pump by turning the lever counterclockwise until it clicks into place.

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**Installation and
operating instructions**

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Connecting the waste container
Connect the end of the suction tube which is not
close to the patient (with the text 'Waste') to the
waste container.

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- Header
- Footer
- Content area

Buttons

① Shortcut to the 'Settings' submenu. Here, amongst other things, the following can be set

- the desired screen language
- the desired volume of the acoustic signals
- the 'Fluid Supply Monitor'

② Switch acoustics signals ON/OFF

Desired volume: The desired signal volume can be set in the 'Settings' menu (see ①)

The state 'switched off' is retained until the user deactivates this function or switches the device off and on again.

③ Fluid Supply Monitor. When this button is pressed, the menu 'Fluid Supply Monitor' (monitoring of the amount of fluid in the irrigation bag) is called up, but only if this function has been activated in the menu 'Extended Settings' (see ①).

④ Call up the Procedure List.

Displays

⑤ Texts describing the current menu level

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10.2 Start screen

After switching on the supply voltage, a KARL STORZ logo with the device name is displayed for approx. 5 s (operating system loading time) followed by a white screen for at least 2 s.

This should allow the user to detect any pixel errors on the display.

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10.3 Operation Selection menu

10.3.1 Function and use

After starting, the operation type to be performed can be selected.


The operation types 'Hysteroscopy', 'Resection' and 'Laparoscopy' are available, and can be activated using the buttons ① or ② or ③.

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10.4 Procedure List

10.4.1 Function/use

After selecting the operation type, the Procedure List with the procedures saved for this kind of intervention is displayed. On the first page, the last four procedures to have been used are shown, with the most recently used one first. The top procedure is shown as 'selected' (orange).

As of page 2, the other procedures are listed according to the procedure number.

If the user wishes to create a new procedure, this can be done using the 'Show Sheaths' button or by copying an existing procedure.

The HYSTEROMAT® E.A.S.I. contains numerous predefined procedures. For hysteroscopy, these begin at 150, for resection at 250 and for laparoscopy at 399. Specific key values are saved in each of these procedures.

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

10.4.2 Buttons

- ① Jump to the previous page of the Procedure List
- ② Jump to the next page of the Procedure List
- ③ Shows the current page number / total number of pages.
- ④ Procedure selection is saved. Following this, the system goes to the menu level 'Tube Set Insertion' or into the operating menu 'Normal Operation GYN' depending on whether the tube has already been inserted or not.
- ⑤ Jump back to Operation Selection (Hysteroscopy / Resection / Laparoscopy)
- ⑥ Name of procedures. When you select a procedure, this is highlighted in orange.
- ⑦ Show Sheaths (only for hysteroscopy and resection operations). If you press this button, you can create a new procedure (starting with the selection of a sheath, thus the name 'Show Sheaths')

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10.5 Sheath List

10.5.1 Function/use

! **WARNING:** The HYSTEROMAT® E.A.S.I. relies on the characteristics of the sheath used for correct pressure measurement. Only use the device with gynecology sheaths from KARL STORZ. KARL STORZ cannot accept any liability for the use of sheaths from other manufacturers.

Creation of one's own procedure via sheath selection.

Select one of the sheaths saved for this type of intervention (hysteroscopy / resection) and confirm this selection.

On the first page, the last four sheaths to have been used are shown, with the most recently used one first. The top sheath is shown as 'selected'.

As of page 2, the other sheaths are listed according to the sheath number.

By pressing the button 'Show Procedures' you return to the Procedure List.

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Origin of the values: The values which are displayed in this menu originate from the procedure which was selected after starting.

Example: If an intra-uterine pressure of 40 mmHg was set in the procedure selected after starting, the value 40 mmHg is shown in the menu 'Pressure/Flow Settings'.

Saving changes permanently: Changes which are made here are only retained for as long as no new sheath is selected. If the settings should be retained for longer, the values must be saved in a new procedure. To do this, the button 'Save as New Procedure' ⑤ must be pressed

Overwriting the values of the operating menu:

Any previously changed settings in the operating menu 'Normal Operation GYN' are overwritten upon confirmation of the settings with the button ⑦.

Example: You have manually set an intra-uterine pressure of 35 mmHg in the menu 'Normal Operation GYN'. If you now set a value of 40 mmHg in the menu 'Pressure/Flow Settings' and save this by pressing the button ⑦, from now on a value of 40 mmHg is shown in the menu 'Normal Operation GYN'. This value applies until the device is switched off.

Pressures of 65 mmHg upwards: If the intra-uterine pressure is set to more than/equal to 65 mmHg, a menu is called up once in which it must be expressly confirmed that pressures greater than 65 mmHg are permitted.

Standard specification values for each new procedure:

Intra-uterine pressure 45 mmHg
Normal flow: 50 ml/min
Blood/Mucus Mode Pressure increase: 0%
Blood/Mucus Mode Flow : 150 ml/min

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10.6.2 Buttons

⑧ Increase the specified setpoint for the intra-uterine flow.

Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) which is to be displayed in the 'Intra-uterine pressure' field of the operating menu 'Normal operation GYN' is specified.

⑨ Increase the specified setpoint for the intra-uterine pressure.

⑩ Decrease the specified setpoint for the normal flow.

'Normal flow' refers to the flow associated with Normal Operation and not the flow seen in the High Flow Mode.

⑪ Increase the specified setpoint for the normal flow.

⑫ Decrease the specified setpoint for the percentage increase in the Blood/Mucus Modus.

⑬ Increase the specified setpoint for the percentage increase in pressure in the Blood/Mucus Modus.

⑭ Decrease the specified setpoint for the flow in the Blood/Mucus Mode

⑮ Increase the specified setpoint for the flow in High Flow Mode

⑥ Return to the previous menu; the changed values are discarded

⑦ Return to the previous menu; the changed values are saved

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10.7 Save as New Procedure

10.7.1 Function/ use

The current settings are saved as a new procedure together with a designation, name and number. The designation and name are restricted to 20 alpha-numeric characters or a total width of 200 pixels.

10.7.2 Buttons

- ① Tab – switches between the input windows (display elements 1, 2 and 3)
- ② Return to the menu level 'Pressure/Flow Settings'); the entries are discarded
- ③ Proceed to the menu level 'Normal Operation GYN'; the entries are saved

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10.8 Tube insertion / animation

10.8.1 Function/use

Animation to aid the assistant with insertion of the tube.

This menu level is exited as soon as it is detected that all of the components necessary for the relevant operation (irrigation function / suction function) have been connected.

Example 1: The Continuous Flow Mode is decisive for the selected procedure 'Cont.-Flow HYS B.I.O.H'. For this mode, both an irrigation and suction tube must be inserted. The menu level is only exited once this has been done.

Example 2: For a Single Flow procedure, only the irrigation tube has to be inserted in order to leave this menu level.

Within a fixed period of 24 hours, operation with a day set is permitted if this has been selected accordingly in the 'Extended Settings'.

This means that the irrigation tube (excluding the patient tube element) can remain in the pump.

After switching on the device, a state change in terms of the position detection of the pressure measuring chamber and irrigation pump is not necessary. The relevant animations cease.

The animation is comprised of the following images (see next page).

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10.9.1 Function/use

Explanation of terms: see pages 26-30

Use: Working level for Normal Operation for therapeutic and diagnostic hysteroscopies for the flow modes → 'Continuous Flow' and → 'Single Flow'.

Setting the desired intra-uterine pressure: Using the '+/-' buttons, the setpoint of the displayed cavity pressure can be changed. This can be done independently of the pump state (→activated/ → deactivated).

Blood/Mucus Mode (4) : When the pumps are both activated and deactivated, it is possible to switch between Normal Operation and → Blood/Mucus Mode within Continuous Flow procedure using the button.

Once the Blood Mucus Mode is set, the values which were set in the Pressure/Flow Settings are displayed.

Using the '+/-' buttons, the percentage pressure increase of the intra-uterine pressure can be changed.

The pump cannot be switched on whilst the Blood Mucus Mode button is pressed.

Bipolar Instrument Mode (3) : The → Bipolar Instrument Mode button can only be actuated when the pumps are activated. If a HF electrode or similar instrument is used, this function should maintain the pressure in the cavity.

If no electrode or an instrument which is too small are operated with this function, the control system causes this Mode to be exited automatically to prevent the patient from being injured.

HYS sheaths and RES sheaths have no Bipolar Instrument Mode. This means that this button is not enabled for these sheaths.

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Footswitch: The special modes 'Bipolar Instrument' and 'Blood/Mucus Mode' can be executed either via the touch screen and/or a connected footswitch ⑤. The following applies for the footswitch:

Grey pedal for Bipolar Instrument, blue pedal for Blood Mucus Mode

① NOTE: The footswitch symbol is only displayed if a footswitch is actually connected.

Idle Mode: Idle mode is another fully automated additional function for reducing liquid consumption when changing instruments or when removing tissue samples from the body. The Idle Mode is activated when the flow exceeds a value of 320 ml/min for more than a second.

To ensure patient safety, on activation of the Idle Mode the Instrument Mode is automatically deactivated.

In Idle Mode, the irrigation pump works continuously at a rate of 100 ml/min; the suction pump is motionless to prevent air being sucked in.

The device pressure is shown in this phase. After being in Idle Mode for at least five seconds, once a certain pressure has been reached (depending on the sheath), the display switches from Idle Mode to normal automatic operation.

The changed parameters are retained even after deactivation of the pumps, the Blood/Mucus Mode and/or the 'Bipolar Instrument' Mode switch to Normal Operation.

After switching the line voltage off and on again or by jumping back to the Procedure List, these are then replaced with the values of the selected procedure.

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The number '50', which appears here in the content area, was specified by the manufacturer for this specific procedure.

The predefined number can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer.

Predefinitions: from number 150 for HYS, from number 250 for RES, number 399 [currently only predefined LAP procedure] for LAP

② By pressing the - button the desired intra-uterine pressure in mmHg is reduced.

③ Bipolar Instr.: By pressing this button, the → Instrument Mode (Cut Mode) is switched on. The button is then shown in orange.

④ Blood/Mucus: By pressing this button, the → Blood/Mucus Mode is switched on. The button is then shown in orange.

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⑥ By pressing this button, the pump is activated (the pump reels remove). When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

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10.10 Extended Display GYN

10.10.1 Function and use

The display mode 'Extended Display' which can be enabled via the menu 'Extended Settings' has the same properties and functions as the Normal Mode.

The only difference is that the flow is additionally displayed. This can be adjusted at any time, except with Single Flow sheaths.

10.10.2 Buttons

Explanation of terms: see page 26-30

① By pressing the + button, the desired intra-uterine pressure in mmHg is increased.

The number "50", which appears here in the content area, was specified by the manufacturer for this specific procedure.

The predefined number can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer (from number 100).

② By pressing the - button, the desired intra-uterine pressure in mmHg is decreased.

③ Bipolar Instr.: By pressing this button, the → 'Bipolar Instrument Mode' is switched on. The button is then shown in orange.

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④ Blood Mucus: By pressing this button, the → Blood/Mucus Mode is switched on. The button is then shown in orange.

⑤ Footswitch symbol, if a footswitch is connected

⑥ By pressing this button, the pump is activated (the pump reels move).

⑦ By pressing the + button, the desired flow in ml/min is increased.

⑧ By pressing the - button, the desired flow in ml/min is decreased.

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10.10.3 Displays

- ⑨ Procedure names are comprised of the following 4 elements:
 - Procedure icon
 - Procedure number
 - Designation 1
 - Designation 2
- ⑩ Designation of the pressure/increase in pressure. During Normal Operation: Intra-uterine pressure. When switching on the Blood/Mucus Mode: Increase in pressure (change in percent)
- ⑪ Symbol for pressure
- ⑫ Pressure (during Normal Operation three digits, in Blood/Mucus Mode two digits with percentage sign)
- ⑬ Designation of flow
- ⑭ Symbol for flow
- ⑮ Flow in ml/min

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Confirm sheath!

To start pump confirm connection of the following sheath:
"4mm (26152B/BO)"

Current settings are as follows:

Operation Mode	Hysteroscopy
Sheath	4mm (26152B/BO)
Flow Mode	Continuous Flow
Pressure	40 mmHG
Normal Flow	50 ml/min
High Flow Mode Press. Inc.	0 %
High Flow Mode Flow	150 ml/min

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10. 11 Sheath and settings confirmation

On pressing the start button, you are asked whether the selected sheath is connected. The set or preset pressure and flow values continue to be displayed.

The user is requested to check the settings and sheath and confirm with 'Yes' or 'No'.

On pressing the 'Yes' button ②, the screen jumps to the operating menu ('Normal Operation' GYN' or 'Extended Display GYN') and the pump is automatically started.

On pressing the 'No' button ①, the screen jumps to the procedures and sheaths ('Procedure List' menu).

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
10.12LAP Mode operating menu

10.12.1 Function/use

Use: Normal working level for the laparoscopic mode. Using the \pm buttons, the setpoints for the suction and irrigation flow parameters can be changed. This can be done independently of the pump state (activated/deactivated).

Each individual pump can be activated and deactivated independently with its own Start/Stop button.

10.12.2 Buttons

- ① Decrease the setpoint for irrigation flow, 100 – 1300 ml/min
- ② Increase the setpoint for irrigation flow, 100 – 1300 ml/min
- ③ Activate/deactivate the irrigation pump status display (locking function) with On/Off
- ④ Decrease the setpoint for suction flow, 100 – 1300 ml/min
- ⑤ Increase the setpoint for suction flow, 100 – 1300 ml/min
- ⑥ Activate/deactivate the suction pump status display (locking function) with On/Off

10.12.3 Displays

- ⑦ Laparoscopy mode
- ⑧ Procedure icon, number, designation and name
- ⑨ Text 'Irrigation flow [ml/min]'
- ⑩ Symbol (Irrigation) flow, is always displayed
- ⑪ Flow (four digits, in ml/min)
- ⑫ Text 'Suction flow [ml/min]'
- ⑬ Symbol (suction) flow, is always displayed
- ⑭ Flow (four digits, in ml/min)

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
10.13 Fluid supply monitoring of the irrigation bag

10.13.1 Function/use

Accessing the menu: This menu is accessed by pressing the button 'Fluid Supply Monitor' (see the orange highlighted irrigation bag at the bottom of the screen)

This menu is only displayed if the 'Fluid Supply Monitor' is activated (enabled) in the menu 'Extended Settings'.

Use: This menu monitors the level of fluid in the irrigation bag. This menu is available regardless of whether the pumps are activated or deactivated. For safety reasons, when the pump is activated the display returns to the previous menu if no button is pressed within 3 seconds. Values which were modified will then be discarded.

The selected irrigation bag size and remaining available fluid volume are displayed. The value for the irrigation bag size can be changed with the +/- buttons (buttons ① and ②). Moreover, the remaining volume value can be modified, if required (buttons ③ and ④), for example if the procedure is started with a partly-used irrigation bag. Using the ⑤ button, the available fluid supply can also be reset to the size of the bag, for example when a new full irrigation bag is connected. All functions can be used regardless of the status of the irrigation and suction pumps (activated/deactivated).

- ⑥ Discard the values and exit the mask
- ⑦ Confirm and save the values and exit the mask

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The modified parameters are retained even if the pumps are deactivated. If the power is switched off and on, however:

- the irrigation bag size is reset to the original value which was set in the menu 'Extended Settings'
- the remaining fluid volume is set to the same value as the irrigation bag size (it is assumed that a new irrigation bag is used for each new intervention).

10.13.2 Buttons

① **NOTE:** The values for the irrigation bag size are overwritten by any modifications made in the section 'Irrigation bag size' in the 'Extended Settings' menu.

① **NOTE:** Please do not forget to save your settings **⑦**.

① ② Decrease/increase irrigation bag size setpoint. Possible setpoints: 1 – 3 – 5 – 10.

③ ④ Decrease/increase the actual value for available fluid supply.
Range: 'Irrigation bag size' up to 0.0

⑤ Reset the available fluid supply to the value of the selected bag size (see input field above) - for example when changing the irrigation bag.

⑥ Return to the previous menu level, changes are discarded

⑦ Return to the previous menu level, changes are saved

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11.1 Basic 'Settings' menu

Accessing the menu

This menu is reached by pressing the 'Settings' button ②.

Explanation of terms: see page xx

11.1.1 Function / use

This menu is the basic menu (main menu) in the menu structure 'Settings'. It is the starting point for the submenus in which settings for the Hysteromat® E.A.S.I. can be made and displayed.

11.1.2 Buttons

① Shortcut to the menu 'Pressure/ Flow Settings'. This is where you enter the initial values for the 'operating menu'.

② Shortcut to the menu 'Extended Settings'. In this menu, you can:

- activate the Fluid Supply Monitor
- set the standard value for the irrigation bag size
- activate the use of → Day sets
- activate the "Extended Display GYN" to allow not only the pressure but also the flow to be set.

③ Shortcut to the 'Device Settings' menu. In this menu, you can:

- set the language
- set the display brightness
- set the volume for the information and alarm signals
- set the volume for button actuation
- set the date format
- proceed to the menu 'Date and Time Settings'

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④ Shortcut to the menu 'Device Configuration'. In this menu, the following information is displayed:

- serial number
- operating hours
- user interface software version
- software process control

⑤ Shortcut to the menu 'System Log'. This is where the information and alarm signals are recorded.

⑥ Shortcut to the menu 'Service Programs'. Access only with password for authorized service personnel.

⑦ Shortcut to the menu 'Procedure Administration'. In this menu, you can:

- save predefined procedures (e.g. 100 Cont. Flow HYS B.I.O.H) under another name and set different values for the newly created procedure (e.g. different pressure, different sheath)
- rename procedures (excluding the procedures predefined on delivery).
- delete procedures (excluding the procedures predefined on deliver).
- overwrite procedures with other values (excluding the procedures predefined on deliver).

⑧ Return to the previous menu level.

⑨ Return to the previous menu level. Saving the settings.

NOTE: When you switch to the menu level 'Operation Selection' the buttons ① and ⑦ are not available.

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11.2 Pressure/Flow Settings in GYN Mode

11.2.1 Function/ use

Input menu for the pressure and flow parameters:

- Specification of setpoint for intra-uterine pressure.
Example: In this menu, a value of 40 mmHg is specified. This value of 40 is then saved and displayed in the operating menu 'GYN Mode'. Here, the value for the current operation can be changed at any time.
- Specification of setpoint for normal flow.
Example: In this menu, a value of 40 ml/min is specified. Irrigation with a flow of 40 ml/min. is then achieved in the operating menu 'GYN Mode'. Here, the value for the current operation can be changed at any time.
- Specification of setpoint for pressure increase in High Flow Mode
Example: If a value of 20% is specified in this menu, a pressure increase of 20% in addition to the set intra-uterine pressure (e.g. 120 mmHg) will be realized in the operating menu 'GYN Mode', but only if the 'High Flow Mode' button is pressed (this is only possible in the flow mode 'Continuous Flow' and not in 'Single Flow'). Instead of 120 mmHg, the user would then work with a value of 144 mmHg ($120 + 20\%$ of $120 = 144$). The percentage increase which is specified in this menu can still be changed at any time for the current operation in the GYN operating menu.
- Specification of setpoint for High Flow Mode flow
Example: If a value of 150 ml/min is specified in this menu, the irrigation flow rate will be 150 ml/min during pump operation, but only if the 'High Flow Mode' button was pressed in the operating menu 'GYN Mode' (this is only possible in the 'Continuous Flow' flow mode)

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11.2.2 Buttons

① Decrease the specified setpoint for the intra-uterine flow.

Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) which is to be displayed in the 'Intra-uterine pressure' field of the operating menu 'GYN Mode' is specified.

In the menu 'GYN Mode', this value can still be changed specifically for the operation at any time.

② Increase the specified setpoint for the intra-uterine flow.

③ Decrease the specified setpoint for the percentage increase in pressure in the Blood/Mucus Mode.

④ Increase the specified setpoint for the percentage increase in pressure in the Blood/Mucus Mode.

⑤ Decrease the specified setpoint for the normal flow.

Here 'normal flow' refers to the flow associated with Normal Operation and not the flow seen in the Blood/Mucus Mode.

⑥ Increase the specified setpoint for the normal flow.

⑦ Decrease the specified setpoint for the flow in the Blood/Mucus Mode.

⑧ Increase the specified setpoint for the flow in the Blood/Mucus Mode.

⑨ Return to the previous menu level; the changed values are discarded

⑩ Return to the previous menu level; the changed values are saved

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11.3.1 Function / use

Input menu for the flow parameters:

- Specification of setpoint for irrigation flow.
Example: In this menu, a value of 500 ml/min is specified. This value of 500 is then saved and displayed in the operating menu 'LAP Mode'. The specified value shown here can be changed specifically for the operation at any time.
- Specification of setpoint for suction flow.
Example: In this menu, a value of 500 ml/min is specified. A suction flow of 500 ml/min is then achieved in the operating menu 'LAP Mode'. The specified value shown here can be changed specifically for the operation at any time.

11.3.2 Buttons

- ① Decrease the specified setpoint for the irrigation flow.
Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu; the desired value (setpoint) which is to be displayed in the 'Irrigation flow' field of the operating menu 'LAP Mode' is specified.
- ② Increase the specified setpoint for the irrigation flow.
- ③ Decrease the specified setpoint for the suction flow.
- ④ Increase the specified setpoint for the suction flow.
- ⑤ Return to the previous menu level; the changed values are discarded
- ⑥ Return to the previous menu level; the changed values are saved

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11.4 'Extended Settings' menu

11.4.1 Function/ use

This is an entry menu for further settings.

NOTE: Do not forget to confirm your entries by pressing button ⑦.

11.4.2 Buttons (orange)

① Here you can call up the menu which serves to monitor the irrigation bag fluid supply.

② ③ Decrease/ increase the initial value of the irrigation bag size.

The irrigation bag size is shown in the 'Fluid Supply Monitor' menu if this function has been selected with the button ①.

NOTE: If the initial values are confirmed by pressing the button ⑥, the values set in the 'Fluid Supply Monitor' menu are overwritten!

④ If you press this button, the menu 'Extended Display GYN' is enabled. For a GYN procedure, this allows you not only to specifically set the pressure but also the flow for the current operation.

⑤ Enable the 'Day Set' function.

⑥ Return to the previous menu level and discard the selected values.

⑦ Return to the previous menu level and save the selected values.

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11.5.1 Function / use

This is the menu for setting

- ① ② the language displayed on-screen,
- ⑦ ⑧ monitor brightness
- ③ ④ the volume of information/ alarm signals

! WARNING: Make sure that the set volume of the loudspeakers allows for perception of the acoustic signals despite normal background noise.

- ⑨ ⑩ volume of button actuation
- ⑤ ⑥ date format

In addition, pressing the ⑪ button allows you to switch to the settings menu for date/time. Any modified values that have not been saved are then discarded.

⑫ Return to the basic 'Settings' menu. The modified values are discarded.

⑬ Return to the basic 'Settings' menu. The modified values are saved.

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11.7 'Procedure Administration' menu

11.7.1 Function / use

① By pressing this button, you can overwrite an existing procedure that you defined yourself. Procedures which were defined on delivery cannot be overwritten. This means that the button does not react when you click it for a predefined procedure.

② By pressing this button, you can delete an existing procedure that you defined yourself. Procedures which were defined on delivery cannot be deleted. This means that the button does not react when you click it for a predefined procedure.

③ Save as New Procedure. If you have changed the values of a predefined procedure and want to have these changed values available permanently, we recommend creating a new procedure. To do so, press this button.

④ Rename Procedure. Self-defined procedures can be renamed.

Example: Defining a new procedure and setting other values yourself

Our example has the following objectives:

- A procedure with the sample name '100 Cont. Flow HYS B.I.O.H' is to be copied
- The following value should be changed:
- Pressure: from 50 mmHg to 40 mmHg
- The new procedure should be saved under the name '160 Cont. Flow HYS B.I.O.H 40mmHg' under the number '160'.

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Method for saving a new procedure

- Switch the device on and wait a few seconds until the 'Procedure List' appears.
- Select an operation type.
- On the touch screen, click on the procedure '150 Cont. Flow HYS B.I.O.H'. This procedure is now highlighted in orange. This means the procedure is now activated.
- Click on the tool symbol in the bottom left of the screen.
- (The 'Settings' menu appears).
- Click on the 'Procedure Administration' submenu.
- Click on 'Save as New Procedure'.
- In the orange highlighted field, enter 'Cont. Flow HYS B.I.O.H 40 mmHg'.
- Press the TAB key (the large button on the far bottom left). The second input field is then highlighted orange.
- Leave this field empty
- Press the TAB key again (the large button on the far bottom left). The third input field is then highlighted orange.
- Enter the identification number for the new procedure here, i.e., '160'.
- Click on the check mark in the far bottom right to save your settings.
- (The updated procedure list appears. The new procedure 'Cont. Flow HYS B.I.O.H 40 mmHg' appears at the top and is highlighted in orange, i.e. it is activated)

Method for changing values

- Click on the tool symbol in the bottom left of the screen
- (The 'Settings' menu appears)
- Click on 'Pressure/Flow Settings'

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- Under the header 'Pressure', click twice on the minus sign to reduce the value from 50 to 40.
- Click on the check mark in the bottom right to save your settings.
- (You are once again in the 'Settings' menu)
- (The procedure '160 Cont. Flow HYS B.I.O.H 40 mmHg' has now been permanently changed)
- (The updated procedure list appears. The new procedure '160 Cont. Flow HYS B.I.O.H 40 mmHg' appears at the top and is highlighted in orange, i.e. it is activated)

Method for renaming the procedure

- Click on the tool symbol on the bottom left
- (The 'Settings' menu appears)
- Click on 'Procedure Administration'
- Click on 'Rename Procedure'
- The flashing character (cursor) is now on the left, next to the '1' of '160'.
- Click on the 'arrow-to-the-right button' in the bottom, right-hand area until the cursor is to the right of '40'
- Enter a space. The procedure is now called '160 Cont. Flow HYS B.I.O.H 40 mmHg'
- Save your entries by clicking on the check mark.

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
Method for deleting the procedure

- (You are now in the Procedure List)
- Click on the tool symbol on the bottom left
- (The 'Settings' menu appears)
- Click on 'Procedure Administration'
- Click on 'Delete Procedure'
- Click on the check mark on the bottom right to delete the procedure permanently

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11.8 'Device Configuration' menu

Function/ use

The following elements are displayed:


- Device serial number
- Operating hours
- User interface software version
- Process control software version

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11.9 'System Log' menu

Function/ use

Presentation of the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/information ID.

The most recent entry is allocated the number 0.

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
11.10 'Service Programs' menu

Function/ use

The service programs are only available to service personnel. Access is password-protected.

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Confirm increase of intrauterine pressure!

You are going to increase Intrauterine pressure over 60 mmHg.

Confirm to increase intrauterine pressure over 60 mmHg

Click X to continue with intrauterine pressure of max 60 mmHg

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
11.11 'Confirm increase of pressure' menu

- If the pressure in the activated or deactivated pump status, i.e., also in the 'Pressure/Flow Settings' is increased for the first time to 65 mmHg, this results in a jump to menu level 90, where a reference to the increase in pressure appears.
- This increased pressure is then confirmed once for the relevant procedure using the confirm button ② before the screen jumps back to the previous menu level.
- With the discard button ①, the pressure remains at 60 mmHg and the screen jumps back to the previous menu level. If the pressure is increased again in excess of 60 mmHg, the confirmation process is repeated.

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11.12 'Safe State' menu

- The Safe State is set for certain causes of risks. In this state, the device can no longer be operated normally.
- However, in the 'Settings' area you can access another examination.
- This menu level is a special case. Switching off the line voltage ends the Safe State.

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12 Service and repair

12.1 Fuse replacement

Turn the device off and remove the power cord from the electrical outlet.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Loosen the line fuse holder with a screwdriver or other suitable tool. Remove the line fuse holder.

WARNING: Only use fuses of the correct rating. Replace defective fuses and screw the line fuse holder back on.

Line fuse: 2 x T1.6 AH 250 V

Reconnect the power cord.

Check that the device can be switched on and that the KARL STORZ start screen is displayed.

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12.2 Cleaning and disinfection

! WARNING: Always pull out power plug before cleaning!

! CAUTION: Avoid allowing liquids to enter the housing at all costs.

Wipe-clean the exterior surfaces of the device with a disposable low-lint cloth moistened with disinfectant. For wipe-down disinfection, we recommend using only surface disinfectants based on quaternary ammonium compounds (no solutions predominantly containing alcohol such as rapid-action disinfectants) and "disposable disinfectant methods" based on ammonium chloride.

! CAUTION: When cleaning the device surface, the area around the pressure sensor ⑥ should be wiped carefully due to its sensitive surface.

! WARNING: Do not reuse disposable tubing sets. Discard after use.

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12.3 Maintenance

Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the reliability of the device and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, we recommend a functional or safety test of the device at least once a year. for detailed instructions please refer to the latest version of the service manual.

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FOR DETAILED INSTRUCTIONS PLEASE REFER TO THE LATEST VERSION OF THE SERVICE MANUAL.

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12.4 Service and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine KARL STORZ parts only.

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12.5 Disposal

This device has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

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WARNING: At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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12.6 Repair program

During the repair period, you will generally be loaned a unit, which must then be returned to KARL STORZ as soon as the repaired unit is returned to you.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str.
78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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12.7 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

When shipping the device, use either the original packaging or equivalent packaging marked with the symbols illustrated here.

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12.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet the applicable laws and standards, and
- the device has been used in accordance with its operating instructions at all times.

12.9 Warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot bear the cost of transportation or the risks associated with shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

NOTE: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str.
78532 Tuttlingen, Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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13.1 Information signals

13.1.1 Configuration

All the information signals created, with the exception of those which display the cause for an inoperable device, (entry 'Safe State' in the 'Reaction' column of the table from page 89) are not self-retaining. This means that the signal is only issued whilst the signal conditions are present (unless something different is defined under 'Reaction').

If the signal conditions are only very brief, there are regulations as to the minimum visual/acoustic signal duration to which the described information signals also have to comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by an alarm signal with a higher priority or a low priority information signal is suppressed so long as there is a signal with a higher priority. If there are several signal conditions of the same priority, only the condition which was recognized first will be signaled. Should signals deviate from this rule, this is noted in the corresponding entry from page 89.

13.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, a minimum display duration of 5s is specified.

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13.1.3 Acoustic signals


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. An entire sequence is played through as a minimum.

The acoustic signal can be switched on or off using a button in the footer of the device display.

The acoustic signal, however, can only be switched off temporarily by the user. After switching on the line voltage, the setting is always 'Audio Alarm On'.

Information signals (general)

The pitch is modulated by +/- 2Hz with a modulation frequency of 1.5 Hz. There are 5 different harmonics. In the case of messages which are related to the Safe State, the sequence is repeated every 15s, but otherwise just played through once.

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
Technical description


Specification for acoustic information signals

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

Availability signal


After completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by +/- 2 Hz with a modulation frequency of 1.5 Hz. There are 5 different harmonics.


Blood/Mucus period

The following periodic signal is issued on starting the Blood/Mucus Mode. This signal is repeated every 15 seconds for the duration of the Blood/Mucus Mode.

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


On ending the Blood/Mucus Mode, the following signal sounds:

Low irrigation fluid supply

If the Fluid Supply Monitor is activated and the irrigation fluid supply is low, a double gong is produced. It lasts for 1.5 seconds.

Key click

Pressing a key on the touch screen is confirmed by a short beep. The volume of this beep can be set separately from the volume of all other information signals in the 'Device Settings' menu and can also be permanently deactivated.

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13.2 Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
20: Dual-pedal footswitch recognized Dual-pedal footswitch was plugged in	3	None
21: Dual-pedal footswitch disconnected Dual-pedal footswitch was disconnected	3	End footswitch special function
22: Shaver system recognized The shaver system is connected and switched on	3	None
23: Shaver system disconnected The shaver system was switched off or disconnected	3	End shaver special function
24: Footswitch/Shaver not supported The current operating context has not assigned a function to the actuated foot- or shaver switch	3	None
25: Blood/Mucus mode period expired The Blood/Mucus Mode ended after the set period of time, the footswitch pedal is still actuated	3	None
26: Low irrigation fluid supply The irrigation fluid supply is coming to an end (< 200 ml)	1	None. Message remains until the irrigation fluid supply is above the limit value or the fluid supply monitor is deactivated.
27: Pressure dome detached Pressure dome was detached	2	Cancel activation (if activated)
28: Irrigation pump unlocked Irrigation pump locking lever was opened	2	Cancel activation (if activated)
29: Suction pump unlocked Suction pump locking lever was opened	2	Cancel activation (if activated) Single Flow: no reaction
2A: Illegal password Incorrect password entered	3	None
2B: Password entry mismatch On entering/confirming a new password, different entries were made	3	None
2C: Pressure dome attached On turning on the line voltage, the pressure dome was attached but the day set was not enabled	2	Menu change not executed, message remains until the pressure dome is removed
2D: Invalid procedure number The procedure number which has been selected for saving has already been used	3	None

Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
40: Low capacitance of back-up battery Voltage of the SRAM back-up battery on powering up < 2.7 V	2	None
41: Back-up battery empty (no data storage) Voltage of the SRAM back-up battery on powering up < 2.3 V	1	None
42: Reset of settings/data The settings/data were reset manually or due to an error (one or more saving areas)	2	None
43: Blood/Mucus mode duration exceeded Current Blood/Mucus Mode duration exceeds the setpoint by more than 5s	2	Safe State
44: Software version conflict On switching on the line voltage incompatible software versions were found by the UserInterface and ProcessControl	1	Safe State
45: Activated ext. switch during power-up On switching on the line voltage a pedal on the footswitch or a button on the shaver handpiece was actuated	1	Switch function is not executed, message remains until the switch element is released
46: Inserted external switch is activated On insertion a pedal on the footswitch or a button on the shaver handpiece was actuated	1	Pedal function is not executed, message remains until the switch element is released
47: Activated touch screen during power-up On switching on the line voltage the touch screen was pressed	1	Safe State
48: Touch screen short circuit Maximum touch screen actuation time of 40s exceeded	1	Safe State
49: Pinch-off valve faulty The pinch-off does not function correctly in all positions	2	Cancel activation
4A: Erroneous drive system (irrigation) On switching on the line voltage a faulty component in the drive system was detected	1	Safe State
4B: Erroneous drive system (suction) On switching on the line voltage a faulty component in the drive system was detected	1	Safe State

Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
50: Pressure outside tolerance The value of the exposed sensor is outside the tolerance of 24.4 mV – 427 mV at the sensor output (approx. 4 -70 mmHg)	1	Cancel activation. Pressure dome animation Renewed activation is only possible once the pressure dome has been removed and replaced
53: Critical pressure value A pressure is measured by the sensor which could result in the latter's destruction	1	Cancel activation. Pressure dome animation Renewed activation is only possible once the pressure dome has been removed and replaced
60: Internal temperature too high Temperature inside the device is greater than 55°C	3	None, message remains until the temperature falls below the limit value or is superseded by message 61
61: Critical internal temperature Temperature inside the device is outside the 10°C - 70°C range	2	Cancel activation, message remains until the temperature falls below the limit value
62: Irrigation motor temperature too high Temperature >70°C	3	None, message remains until the temperature falls below the limit value or is superseded by message 63
63: Irrigation motor temperature critical Temperature >100°C or disconnection of a sensor line	2	Cancel activation, message remains until the temperature falls below the limit value
64: Suction motor temperature too high Temperature >70°C	3	None, message remains until the temperature falls below the limit value or is superseded by message 65
65: Suction motor temperature critical Temperature >100°C or disconnection of a sensor line	2	Cancel activation, message remains until the temperature falls below the limit value

Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
70: Sensor test is due after <4h More than 20 hours have elapsed since the last sensor test	2	None, message remains until the pressure dome is detached
71: Sensor test is overdue The offset value of the pressure sensor is no longer up-to-date (pressure dome must be detached and reattached)	1	Activated: Activation remains Message remains Not activated: Pressure dome animation Activation and deletion of the message is only possible once the pressure dome has been removed and replaced
72: Self-test is due after <4h Since the last self-test over 20 hours ago, the device has been in continuous operation	2	None, message remains until the device is switched off and on again
73: Self-test is overdue Since the last self-test over 24 hours ago, the device has been in continuous operation and a new self-test is now necessary	1	Safe State
80: Supply voltage out of range Supply voltage(s) outside the specification	1	Safe State
81: 24 V critical Supply voltage below 15 V	1	Safe State
90: Unable to download a procedure An attempt was made to download a procedure from the RUI whilst the device was not ready-to-receive.	2	None
91: Invalid procedure download data The data sent from the RUI are faulty	2	None
FD: Internal system error I An unexpected system state was detected	3	None
FE: Internal system error II An unexpected system state was detected	2	Cancel activation
FF: Internal system error III An unexpected system state was detected	1	Safe State

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- Type of protection against electric shock: Protection Class I
- Degree of protection against electric shock: Applied part type BF 
- Type of moisture protection: Drip-proof as per IPX1

13.5 Directive compliance (for 26340020-1)

This medical product belongs to Class IIb.

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

13.6 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
13.7 Blockschaltbild

13.7 Block diagram

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**Spare parts,
recommended accessories**

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**14 Spare parts,
recommended accessories**

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15 Appendix

15.1 Cleaning agents and disinfectants

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Typ

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
- (A) = only for flexible endoscopes of the new generation (built later than 2005)

NOTE: You will also find the current list of approved products on our website (www.karlstorz.com).

NOTE: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets

Manual Preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU* telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handpieces, the instructions in their respective instruction manuals must be followed.

The following instruments produced by KARL STORZ are not suitable for cleaning in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass. When preparing and using the solutions, follow the solution manufacturers' instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility by KARL STORZ:

Ecolab GmbH & Co. OHG	Sekusept Plus	1, 2, 3		
	Sekusept Pulver Classic	1, 2, 3		
Estee Pharma GmbH	Ultrasemil AF	1, 2, 3		
	Ultrasemil Aktiv	1, 2, 3		
	Ultrasemil Classic	1, 2, 3		
Fresenius AG	Alid plus Sporcid FF	1, 2 2		
	Hollia Polska So. zo o Indelba Industria E Comercio LTDA	Potasep Hollia Letaldehyde	1, 2 2	
Jose Colado S.A.	Darodol 4000 I knucko Darodol 9000	2 2		
	Darodol Sharkstohd 2000 Aldisan 2000	2 1, 2		
Lysolform Dr. Rosemann GmbH	Almrol Desolform	1, 2 1, 2, 3		
	Lysolform 3000	1, 2		
Medichem International	MediS MediZyme	2 2		
	Edisonite Super Mucadon-1S	1, 3 2, 3		
Merz Hygiene GmbH	Mucadon-2 v-maktiv Mucocil-T neu	1, 3 1, 2, 3		
	MEIASYS Medizintechnik GmbH (formerly IONZA) Orchemie	IO 50 Orchemie	2 1, 3, A	
Promoagent AB	Whycade Descolon Extra Descolon forte	2 1, 2 2		
	Schumacher, Dr. GmbH	Perfektan Endo Perfektan Neu Perfektan TB	2, 3 1, 2 1, 2, 3, A	
Schülke & Mayr GmbH		Gicasept Gicasept AF forte Gicasept FF Gicasept Instru AF Gicasept Med forte Gicazyme Vivacel V	2, 3 1, 2 1, 2 1, 2, 3, A 1, 3 2	
		Sieris	EnvyCare 2 neodisher LM 2 neodisher mediclean neodisher mediclean forte neodisher medizym neodisher SeptoClean	1, 3 1, 3 1, 1 1, 3 1, A 1, 2
			Weigert, Dr. GmbH & Co.	Alcal Plus Matrix
	Whitely Industries PTY.Ltd			

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination shaft, **not** the housing/handle part.

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Typе

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

Note: You will also find the current list of approved products on our website (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handpieces, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles, the respective manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Hersteller/ Manufacturer Производитель	Handelsname/ Commercial name Торговое наименование	Typ/ Type Тип	
Schulke & Mayr GmbH	Thermosept Alka Clean Thermosept DK Thermosept ED Thermosept ER Thermosept R9N-zym	1 2 2 1 1	
	Steris	Instru-Klenz neodisher DuoClean	1, 3 1, A
	Weigerl, Dr. GmbH & Co.	neodisher FA neodisher FA forte neodisher Oxivario neodisher mediclean neodisher mediclean forte neodisher medklar neodisher medlym neodisher SeptoClean neodisher Septo DA neodisher Septo DN	1 1, 3 1 1 1, 3 1 1 1, 2, A 2 2

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15.2 Electromagnetic compatibility (EMC) information

WARNING: Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The HYSTEROMAT® E.A.S.I. model 26340020-1 corresponds to EN/IEC 60601-1-2:2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The HYSTEROMAT® E.A.S.I. model 26340020-1 is a Group 1 device (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the devices to different electrical circuits.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The HYSTEROMAT® E.A.S.I. model 26340020-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the electronic HYSTEROMAT® E.A.S.I. model 26340020-1. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
MKF 4 E-4396	Yes	4	No	Footswitch
Power Cord	No	2	No	Power Supply
SCB Cord	Yes	1.0 m	No	Connecting cable
PE	No	> 3.0 m	No	Potential equalization

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Table 1

Guidance and manufacturer's declaration – electromagnetic emissions

HYSTEROMAT® E.A.S.I. model **26340020-1** is intended for use in the electromagnetic environment specified below.
The customer or user of the HYSTEROMAT® E.A.S.I. model **26340020-1** should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidelines
RF emissions CISPR 11	Group 1	The HYSTEROMAT® E.A.S.I. model 26340020-1 uses only RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The HYSTEROMAT® E.A.S.I. model 26340020-1 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2
Guidance and manufacturer's declaration – electromagnetic immunity

HYSTEROMAT® E.A.S.I. model 26340020-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the HYSTEROMAT® E.A.S.I. model 26340020-1 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 4

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

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Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the HYSTEROMAT® E.A.S.I. model 26340020-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = [3.5/3]\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range*. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 6
**Recommended separation distances between portable and mobile HF communications equipment and the
HYSTEROMAT® E.A.S.I. model 26340020-1**

The HYSTEROMAT® E.A.S.I. model 26340020-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance <i>d</i> [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:
Company stamp/signature:

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseal.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseal.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 5250-5607, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
R. Joaquim Floriano, 413 – 20º andar –
Itaim Bibi, Sao Paulo SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsättrågränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Italahdenkatu 23a, 00210 Helsinki, Finland
Phone: +35-896824774, Fax: +35-896824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopia France S. A.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece E.P.E.
Sokratous & Kyprou 2, 151 27 Melissia, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsiantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for industrial Endoscopy

OOO KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
Solidere – Beirut Souks
Block M, 3rd Floor
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061
Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoskopy Kasachstan
Bokeykhan Ulica 8/2, VP – 1
010000 Astana, Respublika Kasachstan
Phone: +7 7172 57-52 16 / 57 28 49 / 57 09 34
Fax: +7 7172 43 96 96
E-Mail: info@karlstorz.kz

KARL STORZ Endoscopy – Guif & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O.Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.ksegne@karlstorz-rce.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911
Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cnk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
06-07 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 65325448, Fax: +65 65323832
E-Mail: infoesp@karlstorz.com.sg
E-Mail: servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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