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The screenshot shows a software interface for cine mode. At the top, there is a status bar with a camera icon, '99,974 F', '25,997 s', and a clock showing '03:54:02 PM'. Below this is a 'Cine' section with 'Record' options: 'Fluoro', 'HLF', 'Roadmap', and 'Subtract'. The 'Record Rate' is set to '8 fps', with other options being '4 fps', '12 fps', and '25 fps'. The interface is divided into three main panels: 'Landmarking' (LM) on the left with a '+' button and a vertical slider; 'Mask Registration' in the center with directional arrows and a refresh icon; and a right-hand panel with buttons for 'Roadmap Mask', 'View Subtracted', 'Peak Opacify', and 'Set Mask'. A playback control bar at the bottom features a progress slider, play/pause, stop, and fast forward buttons, along with '1F' and '0S' indicators and 'Run 0 of 0'. A 'Minimize' button is on the bottom left and an 'Exit' button is on the bottom right. A taskbar at the very bottom contains several icons, with the video player icon highlighted. Red boxes and lines highlight specific controls, with white text boxes and red text labels explaining their control options.

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
The screenshot shows a software interface for cine mode. At the top, there is a status bar with a camera icon, '99.974 F', '25,997 s', and '03:54:02 PM'. Below this, the 'Cine' section contains 'Record' options: 'Fluoro', 'HLF', 'Roadmap', and 'Subtract'. 'Record Rate' options are '4 fps', '8 fps', '12 fps', and '25 fps'. The 'Landmarking' section has a vertical slider and a '+' button. The 'Mask Registration' section has directional arrows and a refresh icon. On the right, there are buttons for 'Roadmap Mask', 'View Subtracted', 'Peak Opacity', and 'Set Mask'. A playback control bar at the bottom features a progress slider, play/pause, stop, and next buttons, along with '1F' and '0S' indicators and 'Run 0 of 0'. A 'Minimize' button is on the left and an 'Exit' button is on the right. A taskbar at the very bottom contains several icons, with the video camera icon highlighted.

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
The screenshot displays the Cine module interface with the following elements and annotations:

- Top Bar:** Shows system information including 99.974 F, 25.997 s, and 03:54:02 PM. An annotation points to this area: "Control by OptionCineX(X=8 or 12 or 25)".
- Record Section:** Includes checkboxes for Fluoro, HLF, Roadmap, and Subtract. An annotation points to the Roadmap checkbox: "Control by OptionDSARoadmap".
- Record Rate Section:** Features radio buttons for 4 fps, 8 fps, 12 fps, and 25 fps. An annotation points to the 12 fps and 25 fps options: "Control by OptionCine25".
- Landmarking Section:** Contains a vertical slider and a plus sign button. An annotation points to this section: "Control by OptionDSARoadmap".
- Mask Registration Section:** Includes directional arrow buttons and a refresh button. An annotation points to this section: "Control by OptionCine12".
- Right Panel:** Contains buttons for Roadmap Mask, View Subtracted, Peak Opacity, and Set Mask. An annotation points to the Peak Opacity button: "Control by OptionCineX(X=8 or 12 or 25)".
- Timeline Section:** Features a progress bar, playback controls (stop, play, next), and a "Run 0 of 0" indicator. An annotation points to this section: "Control by OptionCineX(X=8 or 12 or 25)".
- Bottom Bar:** Includes Minimize and Exit buttons, and a row of icons. An annotation points to the video camera icon: "Control by OptionCineX(X=8 or 12 or 25)".

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COPYRIGHT © 2019 BY GE HEALTHCARE
RX ONLY - U.S.A
GE HUALUN MEDICAL SYSTEMS CO. Ltd