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**INSTRUCTION**  
for use of medical devices

“Oscillating saw blade, sterile, disposable for Endoprosthetics” (see Appendix)

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Managing Director

1. **Name of medical device**

Oscillating saw blades, sterile, disposable for Endoprosthetics (hereinafter referred to as "saw blade") are designed in accordance with modern requirements to safety and quality and conform to the following standards:

Certificates:

ISO 9001:2015

EN ISO 13485:2016

EU Certificate Directive 93/42/EEC Appendix II, excluding Section 4

Manufacturer: Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25, 32657 Lemgo, Germany.

Class of potential risk of medical device application: 2a

The instruments are made of high-quality stainless steel 1.4028Mo and 1.4031 Mo.

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Fig.12

Fig.13

- Processing of foot / knee / hip joints during endoprosthesis replacement

List of diseases requiring treatment with the use of oscillating saw blades according to clinical judgment:

- **Defective positions, deformations**, for example, deformation of hip joint, asymmetry / flat foot / dislocation / valgus legs, hallux valgus / bent limbs inside / disease of the first metatarsophalangeal articulations.
- **Fractures**, for example, humerus fracture
- **Diseases**, for example, rheumatoid and posttraumatic arthrosis (coxarthrosis, gonarthrosis), arthrosis (osteoarthritis, coxarthrosis, gonarthrosis), rare diabetic arthropathia of leg, Charcot foot, necrosis of ankle bone, arthragra, metatarsalgia

**A surgeon or surgical assistant is responsible for selection of the corresponding instruments.**

### 3. **Contraindications. Adverse effects.**

All general contraindications for orthopedic, surgical interference shall be taken into consideration.

- 1) Metal allergy and hyper sensibility to foreign substances
- 2) Insufficient bone mass and bony tissue quality
- 3) Acute infection
- 4) Operations in which the use of a motorized saw involves risks outweighing the potential benefits. A motorized saw should not be used if any other contraindications apply that prohibit orthopedic surgery in general. Contraindications of the implant must be respected.

Examples of contraindications are not associated directly with proper and correct use of the corresponding medical devices. Proper allowance must be made here for different factors that may be revealed and put in perspective only by a responsible doctor or operator.

- Do not use the instruments for purposes other than that intended.
- It is forbidden to work with the instruments that do not meet the requirements that are imposed thereon.
- The instruments with the following defects shall be immediately rejected: dulled or broken saw blade; deformed (for example, curved or twisted instruments).

**A surgeon or surgical assistant is responsible for selection of the corresponding instruments.**

**When using the instruments, adverse effects are not found.**

### 4. **How to Use.**

The saw blades are meant for single use and supplied in a sterile package. Dispose after use in accordance with the existing regulations.

Instruments may only be used in the appropriate handpiece/power systems with connection types: DePuy Synthes/Zimmer Biomet, Stryker, Aesculap 3Ti. The instructions for use for the power system used must be observed. Use only those saw blades that are adapted to your equipment.

The instrument should be carefully inspected before use to ensure full functionality. Make certain that the saw blades are sterile and ideal regarding their technical and sanitary condition - the sterile package is not impaired; the instruments are not damaged. Consequently, the instruments need to be discarded if the expiry date indicated on the product label has been exceeded, the packaging is opened or damaged or in case of illegible markings. The saw blades shall be fixed so much deep in the tip of the surgical saw as possible. The saw blades that are not inserted to their full length of clamp account for appearance of significant centrifugal force and vibration. Prior to use, make sure that the instrument is safely inserted and locked in its proper position. Contact of the saw blade with metal and other metallic objects must be avoided. The saw blades should be actuated in/at the saw guide and before contact is made with the bone. The thickness of the saw blade must be inferior to that of the slot in the sawing template. Before use on the patient, free mobility of the saw blade in/at the template must be assured. It is necessary to avoid jamming of the saw blade, using it for leverage. This leads to excessive heating, jamming of the saw blade or thermal necrosis as well as to the risk of fracture. It is necessary to avoid a contact of the saw blade with a template, guiding device, or other metal objects while performing oscillating movements. This leads to damage to the saw blade or tools as well as blunting of the blades, injury to bones and surrounding tissue (incl. nerves and blood vessels). This also increases the risk of inadvertent creation of metal chips. In case of breakage of metallic instruments fragments and/or components can be located by medical imaging procedures. Improper use leads to out-of-specification results and injury risk.

Excessive contact pressure must be avoided as this may lead to:

- damage to the tooth surfaces
- reduced service life of saw blade
- excessive heat
- thermal necrosis
- fracture of the blade
- jamming of the blade inside the bone
- undesirable rough and uneven surface

It is necessary to provide sufficient cooling (saline solution) during use. Insufficient cooling leads to clogging of the teeth with bone chips, thermal necrosis, reduction of the service life. Check saw blades for damages and dulling of teeth by means of magnifying arrangements.

Examine the instruments specifically for:

- splits and dulling (rounding) of teeth
- damage in the area of fixation
- curvature of a saw blade

Damaged saw blades (scratches, dents, blunt blades, bent instruments) may lead to the fracture of the instrument or treatment failure such as the patient getting injured by metal fragments left behind in the operative area.

**Nonobservance of these recommendations may lead to the instrument breakage and pose an increased risk for a patient and surgeon.**

## 5. Measures of Precaution in Using.

The saw blades are meant for single use and supplied in a sterile package. Dispose after use in accordance with the existing regulations.

Do not use the instruments if the expiry date indicated on the product label has been exceeded, the packaging is opened or damaged, in case of illegible markings or corrosion. Do not use the instrument if damaged or in case of blunt cutting edges.

Instruments may only be used in the appropriate handpiece/power systems with connection types: DePuy Synthes/Zimmer Biomet, Stryker, Aesculap 3Ti. The instructions for use for the power system used must be observed. Use only those saw blades that are adapted to your equipment.

The saw blades shall be fixed so much deep in the tip of the surgical saw as possible. The saw blades that are not inserted to their full length of clamp account for appearance of significant centrifugal force and vibration. Prior to use, make sure that the instrument is safely inserted and locked in its proper position. Contact of the saw blade with metal and other metallic objects must be avoided. The saw blades should be actuated in/at the saw guide and before contact is made with the bone. The thickness of the saw blade must be inferior to that of the slot in the sawing template. Before use on the patient, free mobility of the saw blade in/at the template must be assured. It is necessary to avoid jamming of the saw blade, using it for leverage. This leads to excessive heating, jamming of the saw blade or thermal necrosis as well as to the risk of fracture. It is necessary to avoid a contact of the saw blade with a template, guiding device, or other metal objects while performing oscillating movements. This leads to damage to the saw blade or tools as well as blunting of the blades, injury to bones and surrounding tissue (incl. nerves and blood vessels). This also increases the risk of inadvertent creation of metal chips. Improper use leads to out-of-specification results and injury risk.

Excessive contact pressure must be avoided as this may lead to:

- damage to the tooth surfaces
- reduced service life of saw blade
- excessive heat
- thermal necrosis
- fracture of the blade
- jamming of the blade inside the bone
- undesirable rough and uneven surface

It is necessary to provide sufficient cooling (saline solution) during use. Insufficient cooling leads to clogging of the teeth with bone chips, thermal necrosis, reduction of the service life. Shifting of the template or slipping of the saw blade can lead to cutting mistakes. Check saw blades for damages and dulling of teeth by means of magnifying arrangements.

Damaged saw blades (scratches, dents, blunt blades, bent instruments) may lead to the fracture of the instrument or treatment failure such as the patient getting injured by metal fragments left behind in the operative area.

Due to wear of the material, a safe and risk-free reuse of the product cannot be guaranteed. Reprocessing of single-use device may create a risk of contamination and can result in infection and irreversible injury of the patient.

The instruments are not provided with a measuring function or an external or internal energy source and do not incorporate software. The instruments do not contain any medicinal product constituents,

ingredients of animal origin or substances or derivatives made of human blood. The devices do not emit (ionizing) radiation or create magnetic fields. As regards to this, no risks are to be expected.

**Nonobservance of these recommendations may lead to the instrument breakage and pose an increased risk for a patient and surgeon.**

## 6. Product specification.

This medical device belongs to the group of devices for professional use.

### Shape of a working part:

#### Parallel.

Dimensions:

Width of a working part: 13.0; 19.05; 25 mm.

Thickness of a working part: 0.89; 1.0; 1.19; 1.25; 1.37; 1.47 mm.

Length of a working part: 60.0; 90.0 mm.

#### Flared.

Dimensions:

Width of a working part: 13.0; 21.0 mm.

Thickness of a working part: 1.19; 1.25; 1.47 mm.

Length of a working part: 70.0; 90.0 mm.

- Surgical motor systems, which the instruments are used with for all common surgical motor systems with connection types:

- Aesculap 3Ti
- DePuy Synthes/ Zimmer Biomet
- Stryker

For example:

- Trauma Recon System, Battery Power Line II, manufacturer DePuy Synthes, RC ФС3 2012/12800
- Zimmer Universal Power System, manufacturer Zimmer Surgical SA, RC P3H 2013/254
- Sodem with an oscillating hand piece Zimmer, manufacturer Sodem Diffusion S.A. RC ФС3 2007/00643
- Stryker System G, manufacturer Stryker Instruments. RC P3H 2017/5929
- Acculan 3Ti, manufacturer Aesculap AG. RC ФС3 2007/00666, ФС3 2007/00665

- Connection with a surgical motor system: inserted into the fastening node of the surgical motor system tip up to the maximum depth until it stops

- Dulling (wearing/grinding): single use

- Shifting of teeth – no shifting

Table 1. Specification of a working part of saw blades

Article	Dimensions (+/- mm)										Shape of working part	Material	Material grade
	Total length (mm)	Tolerance for Total length (mm)	Length of working part (mm)	Tolerance for length of working part (mm)	Width of cutting working part (mm)	Tolerance for width of cutting working part (mm)	Thickness of cutting edge of the working part (mm)	Tolerance for thickness of cutting edge of the working part (mm)	Thickness of saw blade, not including shifted tooth structure (mm)	Tolerance for thickness of saw blade, not including shifted tooth structure (mm)			
05.003.400S	84.20	+/-0,3	60.00	+/-0,3	13.00	+/-0,2	0.89	+/-0,02	0.89	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.401S	114.20	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	0.89	+/-0,02	0.89	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.402S	114.20	+/-0,3	90.00	+/-0,3	19.05	+/-0,2	0.89	+/-0,02	0.89	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.403S	114.20	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	0.89	+/-0,02	0.89	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.411S	114.20	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.00	+/-0,02	1.00	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.413S	114.20	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.00	+/-0,02	1.00	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.421S	114.20	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.422S	114.20	+/-0,3	90.00	+/-0,3	19.05	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.423S	114.20	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.424S	114.20	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.425S	114.20	+/-0,3	90.00	+/-0,3	21.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.431S	114.20	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.432S	114.20	+/-0,3	90.00	+/-0,3	19.05	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.433S	114.20	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.434S	94.20	+/-0,3	70.00	+/-0,3	13.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.435S	114.20	+/-0,3	90.00	+/-0,3	21.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.442S	114.20	+/-0,3	90.00	+/-0,3	19.05	+/-0,2	1.37	+/-0,02	1.37	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo

05.003.443S	114.20	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.37	+/-0,02	1.37	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.451S	114.20	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.452S	114.20	+/-0,3	90.00	+/-0,3	19.05	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.453S	114.20	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.454S	94.20	+/-0,3	70.00	+/-0,3	13.00	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.500S	92.10	+/-0,3	60.00	+/-0,3	13.00	+/-0,2	0.89	+/-0,02	0.89	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.521S	122.20	+/-0,5	90.00	+/-0,5	13.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.523S	122.15	+/-0,5	90.00	+/-0,5	25.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.524S	122.20	+/-0,5	90.00	+/-0,5	13.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.525S	122.20	+/-0,5	90.00	+/-0,5	21.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.533S	122.15	+/-0,5	90.00	+/-0,5	25.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.534S	102.15	+/-0,3	70.00	+/-0,3	13.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.535S	122.15	+/-0,5	90.00	+/-0,5	21.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.553S	122.15	+/-0,5	90.00	+/-0,5	25.00	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.554S	102.15	+/-0,3	70.00	+/-0,3	13.00	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.703S	117.70	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	0.89	+/-0,02	0.89	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.721S	117.70	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.723S	117.70	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.19	+/-0,02	1.19	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.733S	117.70	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo
05.003.734S	117.70	+/-0,3	90.00	+/-0,3	13.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.735S	117.70	+/-0,3	90.00	+/-0,3	21.00	+/-0,2	1.25	+/-0,02	1.25	+/-0,02	Flared	Stainless steel	1.4031Mo/1.4028Mo
05.003.753S	117.80	+/-0,3	90.00	+/-0,3	25.00	+/-0,2	1.47	+/-0,02	1.47	+/-0,02	Parallel	Stainless steel	1.4031Mo/1.4028Mo

Material of a working part and shank: stainless steel 1.4028Mo and 1.4031Mo.

Manufacturer of stainless steel 1.4028Mo:

Zapp Precision Metals GmbH, Letmather Straße 69, 58239 Schwerte, Germany

Manufacturer of stainless steel 1.4031Mo:

voestalpine Precision Strip GmbH, P.O. Box 503, 68428 Munkfors, Sweden  
and Sandvik Materials Technology EMEA AB, Box 12, 16493 Kista, Sweden

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[www.goszdramnadzor.gov.ru](http://www.goszdramnadzor.gov.ru)

Table 2. Specification of saw blade shanks

Shank type	Number of shank in the article	Specimen of instrument article	Shank material	Roughness (Ra)	Hardness (HV)	Shank view
Aesculap 3Ti	7	05.003.723S	Stainless steel 1.4031Mo/1.4028Mo	0.315 - 0.532	510 - 580	 Aesculap® Acculan 3Ti
DePuy Synthes / Zimmer Biomet	4	05.003.422S	Stainless steel 1.4031Mo/1.4028Mo	0.315 - 0.532	510 - 580	 DePuy Synthes / Zimmer®
Stryker	5	05.003.521S	Stainless steel 1.4031Mo/1.4028Mo	0.315 - 0.532	510 - 580	 Stryker®

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## 7. Requirements to Environmental Conservation

The device is not toxic for the environment.

### **Environmental risks of medical devices:**

Corresponding influences are not present.

## 8. Disinfection, cleaning and sterilization

The saw blades are meant for single use and supplied in a sterile package. Dispose after use in accordance with the existing regulations.

## 9. Sterility

The devices are sterilized by means of gamma-radiation **VDmax<sup>25</sup>**. The sterilization is performed in accordance with the requirements of ISO 11137-2:2013.

There is no influence of manufacturing methods (including sterilization) on materials characterization and working properties. The required sterilization method does not jeopardize its safety and operational characteristics.

The sterilization process is validated and subject to process inspection as evidenced by the Sterilization Validation Report.

The devices supplied as sterile are indicated with a mark "STERILE" .

The devices shall be removed from the package under aseptic conditions. The manufacturer may not ensure sterility if the package integrity is damaged, or the package is opened incorrectly, and is not responsible in the circumstances. In this case, the instruments shall be disposed in accordance with internal laws of the Russian Federation (SanPiN 2.1.7.2790-10). Hazard class of medical wastes B.

## 10. Applied International Standards.

The devices "Oscillating saw blades, sterile, disposable for Endoprosthetics" conform to the requirements of the following standards:

DIN EN ISO 13485

DIN EN ISO 10993

DIN EN ISO 11137

DIN EN ISO 14971

## 11. Notification of unwanted events/occurrences:

All serious incidents occurring in connection with the product shall be reported to the manufacturer and to the competent national authority.

## 12. Risk Management

The purpose of risk management is to determine if the risk management process can eliminate or reduce risks for the manufacturer as far as possible, to take measures to reduce risks which cannot be eliminated, and to inform users about residual risks.

The risk management system of Gebr. Brasseler GmbH & Co. KG is intended for identification of hazards that are related to the device / manufacturing process and for assessment of risk related thereto, for control of these risks and keeping track of efficiency of control during the whole service life of the device. Risk analysis, risk assessment begins in the process of design and development. The efficiency of every device in relation to these and any unexpected risks is controlled as part of the management system of Gebr. Brasseler GmbH & Co. KG and its program of post-manufacturing control.

Application of the risk management process meets the requirements of the international standard ISO 14971-2007.

### 13. Storage and Application Regulations

The instruments are packed under sterile conditions.

#### Storage conditions:

Prior to use, the instruments shall be stored in the original packing, away from direct sunlight, dust, moisture, sources of recontamination and high temperature. Do not store the instruments in the same room with reagents and chemicals.

All works with contaminated instruments shall be performed only in safety gloves.

#### Operating Conditions:

Temperature: from +5°C to +45°C

Relative air humidity: from 30% to 80%

The instruments shall be stored in dry condition at the temperature from -50°C to +50°C and at humidity from 10 to 95%.

The instruments with an impaired sterile package as well as damaged instruments shall be rejected immediately.

### 14. Marking of the Medical Device

The following marking is used in the medical device:

Table 3. Marking of the Medical Device.

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## 15. Supply Package

The instruments are supplied in a sterile package, 1 piece per package, with an instruction for use inside of the package.

## 16. Package of the Medical Device

The package is sterilized by means of gamma-radiation **VDmax25**.

The saw blades are packed in a sterile package, 1 piece per package. The package consists of the external wrapping material **GC1** (cardboard "Incada Silk", 260 g / m<sup>2</sup>), that includes saw blades packed in the blister card (package dimensions: 103 x204 mm), and of the sterile double barrier system (internal blister card 170 (±0,5 mm) x 50 (±0,3 mm) mm and external blister card 200 (±0,5 mm) x 95 (±0,3 mm) mm).

Material, of which the internal blister card is made, is a clear film **PA/PEM 70/90** and sealing film **Ecopet P 400**.

Material, of which the external blister card is made, is a clear blister film **PETG / Pentamed® MD-EG00R01-51G8850-100\_0-ZZZ** and sealing film **Tyvek 1073B / HCW CR27**.

The paper label is a white writing paper.

The sterilization validity makes up maximum 5 years from the date of manufacturing specified on the package. The instruments shall be disposed if the minimum shelf life is exceeded, the package is opened or damaged, or data are illegible. The patient labels must be added to the patient's health record, hospital record etc. after initial opening.

Each package is accompanied by a leaflet (Instructions for use).

Technical Specifications of PA/PEM 70/90

Properties	Values	Measurement units	Method
Thickness	160 (± 10 %)	µm	DIN 53370 (arithmetic mean is calculated)
Density	162 (± 8 %)	g/m <sup>2</sup>	ISO 2286-2
Recommended sealing temperature	> 125	°C	Test carried out by the supplier of Südpack
Oxygen permeability	<= 11	cm <sup>3</sup> /m <sup>2</sup> per day bar	23°C, relative air humidity 50 % DIN 53380
Water vapor permeability	<= 4	g/m <sup>2</sup> per day	38°C, relative air humidity 90% ASTM F 1249
Width	See the order	mm	SPO
Roll length	See the order	m	SPO

Technical Specifications of Ecopet P 400

Properties	Values	Measurement units	Method
Thickness	400 (± 5 %)	µm	DIN 53370 (arithmetic mean is calculated)
Density	516 (± 8 %)	g/m <sup>2</sup>	ISO 2286-2
Recommended sealing temperature	> 120	°C	Test carried out by the supplier of Südpack
Oxygen permeability	<= 10	cm <sup>3</sup> /m <sup>2</sup> per day bar	23°C, relative air humidity 50 % DIN 53380
Water vapor permeability	<= 5	g/m <sup>2</sup> per day	38°C, relative air humidity 90% ASTM F 1249
Width	See the order	mm	SPO
Roll length	See the order	m	SPO

ASTM F1249 = standard procedure of velocity determination of water vapor transmission through a plastic film and film with using a modulated infrared sensor.

SPO = single point displacement

Technical Specifications of PETG / Pentamed® MD-EG00R01-51G8850-100\_0-ZZZ

Properties	Standard	Measurement units	Values
Thickness	DIN 53370 / ISO 4593	µm	200 ... 1200
Density	DIN EN ISO 1183-2	g/cm <sup>3</sup>	1.26 ... 1.30
Tensile strength, min.	DIN EN ISO 527 speed of testing V50 mm/min, measured along the length, depending on thickness	MPa	nevertheless 40
Dynamic strength, min.	DIN EN ISO 8256 measured along the length	kJ/m <sup>2</sup>	nevertheless 300
VICAT-melting point	DIN EN ISO 306 measured in oil, method B/50	°C	78 ... 82
Tolerance for width	DIN 53570	mm	-1 ... 1


Technical Specifications of Tyvek 1073B / HCW CR27

Properties	Method	Value
Basis weight	ISO 536	75 (± 4 %) g <sup>-2</sup>
Weight of coating	ASTM F2217:02	13.2 g <sup>-2</sup> Min. 12.6 g <sup>-2</sup> Max. 14.8 g <sup>-2</sup>
Gross weight	ISO 536	88.2 (± 4 %) g <sup>-2</sup>
Tensile strength (MD <sup>1</sup> ) longitudinal	ISO 1924-2	7.6 N/mm Min.: 6 Max.: 9.2
Tensile strength (CD <sup>2</sup> ) transverse	ISO 1924-2	8.2 N/mm Min.: 6 N/mm Max.: 10.4 N/mm
Ultimate elongation (MD <sup>1</sup> ) longitudinal	ISO 1974	3425 MN Min.: 2135 MN Max.: 5160 MN
Ultimate elongation (CD <sup>2</sup> ) transverse	ISO 1974	3514 MN Min.: 2313 MN Max.: 5160 MN
Air permeability	ISO 5636/3	140 ml.min. <sup>-1</sup> Min.: 90 ml.min. <sup>-1</sup> Max.: 560 ml.min. <sup>-1</sup>
Strength of seam	ASTM F88	nevertheless 4N /25mm
Width of seam	ISO 11607-1	5 (0/+ 1 mm) mm

The following information is specified on the package of instruments:

- logo and name of the company, for which the instruments are manufactured

- name of trademark which applies for all saw blades DePuy Synthes, intended for conducting operations on reconstruction of joints.

- catalogue number

- lot number

**Oscillating Saw Blade - for DePuy Synthes® / Zimmer®**

- name of the instrument and connection type of the surgical system, which the instrument is used with

**90x13x1.19mm**

- dimensions of the instrument: 90 - length of working part, 13 - width of cutting working part, 1,19 - thickness of saw blade, not including shifted tooth structure

- GS1 bar-code


- number of pieces in the package

- do not use if the package is damaged

- The symbol points out the necessity for the user to familiarize himself with important information in the instruction for use, such as warnings and precautionary measures which for a variety of reasons may not be placed on the medical device.

2024-09-31 - date after which the medical device may not be used any more

2019-04-03 - date of manufacture

- indication that this device is sterilized by means of radiation

 - only for single use

 - symbol "Refer to the operating manual"

**Rx Only** - indication that the US Federal Laws limit the sale of this device by a doctor or at the request of a doctor

- GS1 Data matrix code

The symbols glossary is provided within the instruction for use

- description of symbols is provided in the instruction for use

**MAT** Stainless Steel

- material of which the instrument is made is stainless steel

(01) 14053613217568  
(11) 190403  
(17) 240331  
(10) 819080

- This is a plain text of GS1-code. The code points out our basic primary UDI = unique device identification.

BG: нож за осцилираща глава	LT: Vibruojanti pjūklis galvė
HR: Oscilirajući list pile	NO: Oscillerende sagblad
CS: Oscilující pilový list	PT: Lâmina de serra oscilante
DA: Oscillerende savklinge	RO: Lamă de ferstrău oscilantă
EN: Oscillating Saw Blade	RU: Осциллирующее пильное полотно
FR: Lame de scie oscillante	SR: Oscilirajući list testere
DE: Oscillierendes Sägebblatt	SK: Oscilujúci pilový list
EL: Πολυτρομική Λαμίδα πριονιού	SL: Oscilirajoč žagin list
HU: Oscilláló fűrészlap	ES: Hoja de sierra oscilante
IT: Lama per sega oscillante	SV: Oscillerande sågblad
LV: Oscilējošais zāģa asmens	TR: Salınimli testere bıçağı
NL: Oscillerend Zaagblad	PL: Oscylujący brzoźczoł pily

- name of the instrument in different languages

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- information about the manufacturer

- layout view of the instrument

- indication of CE certificate availability

**DIST**  
Synthes GmbH  
Emattstrasse 3  
4436 Oberdorf  
Switzerland

- name of distributor

Distributed in US by:  
The Anapach Effort, Inc.  
4590 Riverside Drive  
Palm Beach Gardens, FL 33410  
USA

- name of distributor in the USA

Made in Germany

- country of origin

On the blister, in which the instrument is placed, the following information is specified:

**REF** - catalogue number

**LOT** - lot number

**Oscillating Saw Blade** - name of the instrument

**2019-04-03** - date of manufacture


**2024-03-31** - date after which the medical device may not be used any more

**GS1**  
- GS1 Data matrix code

**(01) 14053813217868**  
**(11) 190403**  
**(17) 240331**  
**(10) 819080** - This is a plain text of GS1-code. The code points out our basic primary UDI = unique device identification.

**DePuy Synthes** - logo and name of the company, for which the instruments are manufactured

The following information is specified on the Russian-language label of the package of instruments:

**REF** - catalogue number

Label ver.0.20 - QR code of the label

- only for single use

- indication that this device is sterilized by means of radiation

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
Width of a working part (mm)	Thickness of cutting edge of the working part (mm)	Connection type	Designation of connection type
	0,89 mm - 09		7 – Aesculap 3Ti
	1,00 mm - 10		4 – DePuy Synthes/ Zimmer Biomet
	1,19 mm - 12		5 – Stryker
	1,27 mm - 13		
	1,37 mm - 14		
	1,47 mm - 15		

05.003.400S  
 05.003.411S  
 05.003.422S  
 05.003.435S  
 05.003.443S  
 05.003.454S  
 05.003 – series

4 – Connection type: DePuy Synthes/ Zimmer Biomet

0; 1; 2; 3; 4; 5 – thickness of cutting part – 0.89; 1.00; 1.19; 1.27; 1.37; 1.47 mm

0; 1; 2; 3; 4; 5 – width and length of working part, (see definition on the layout above)

S – sterile

05.003.500S

05.003.521S

05.003.533S

05.003.535S

05.003.554S

05.003 – series

5 – Connection type: Stryker

0; 2; 3; 5 – thickness of cutting part – 0.89; 1.19; 1.27; 1.47 mm

0; 1; 3; 4; 5 – width and length of working part, (see definition on the layout above)

S – sterile

05.003.703S

05.003.721S

05.003.734S

05.003.735S

05.003.753S


05.003 – series

7 – Connection type: Aesculap 3Ti

0; 2; 3; 5 – thickness of cutting part – 0.89; 1.19; 1.27; 1.47 mm

1; 3; 4; 5 – width and length of working part, (see definition on the layout above)

S – sterile

1. The width of a working part of the saw blade is measured along the outside edge of two outermost teeth
2. Maximum depth, to which the saw blade can get through the bone
3. The thickness of a cutting working part is measured along the cutting edge
4. Thickness of saw blade, not including shifted tooth structure

Transport package:

A package used for transportation, is chosen specific to every order. A corrugated cardboard is usually used. The dimension of a transport package depends on the volume of goods.

The following dimensions of cardboard boxes are available currently.

180 x 103 x 53
190 x 120 x 80
260 x 140 x 90
280 x 190 x 130
315 x 230 x 125
370 x 280 x 230
600 x 400 x 250
700 x 450 x 300
620 x 620 x 460
350 x 300 x 290
300 x 185 x 290

Marking is applied only to the package of instruments itself. There is no special marking on the transport package. After the devices have been packed at the factory, the transport package is closed carefully and provided with a cargo hauling label and transportation documents. Then the packed devices are transferred to a forwarder and shipped to a customer.

#### 17. Transportation Conditions

The instruments “Oscillating saw blade, sterile, disposable for Endoprosthesis” are transported by all transport modes in covered vehicles according to the rules of carriage applying to this transport mode at the temperature from  $-50^{\circ}\text{C}$  to  $+50^{\circ}\text{C}$  and humidity from 10 to 95%.

The instruments shall not come into direct contact with fire and water, especially with respect to their package. The packed sterile instruments shall be protected from dust, moisture and contamination during transportation and storage. There are no other special conditions of storage, transportation or use that shall be considered.

#### 18. Service life

The sterilization validity makes up maximum 5 years from the date of manufacturing specified on the package

*Control and functional check:*

The instruments having the following defects shall be immediately rejected:

- Deformations (for example, bended/intorted instruments)
- Corrosive surface

*Anticipated Usage Frequency of Saw Blades:*

The saw blades are meant for single use:

- one operation on one patient

Single-use devices (labelled with a sign  ) shall not be reused!

The reuse of such devices leads to the risk of infecting, and safety of devices may not be ensured any more.

## 19. Disposal and Destruction

Disposal of failed instruments is performed in accordance with internal laws of the Russian Federation (SanPiN 2.1.7.2790-10). Hazard class of medical wastes B.

Potentially contaminated, pointed and sharp items shall be disposed in single-use containers.

## 20. Manufacturer's Warranties

The user is responsible for checking the product prior to use to ensure that it is suitable for the intended purpose. The user is responsible for the use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

According to general provisions and sales conditions, the warranty claims of the customer being a commercial entrepreneur become effective on the date of goods delivery and extend for one (1) year. The guaranteed shelf life - sterilization validity makes up maximum 5 years from the date of manufacturing specified on the package.

On issues relating to quality, address "MedTestConsulting" LLC.

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Managing Director

Appendix 1 – List of items

Appendix 2 – Authorized Representative

**“Oscillating saw blade, sterile, disposable for Endoprosthetics”**

**I. Versions:**

05.003.400S; 05.003.401S; 05.003.402S; 05.003.403S; 05.003.411S; 05.003.413S; 05.003.421S;  
05.003.422S; 05.003.423S; 05.003.424S; 05.003.425S; 05.003.431S; 05.003.432S; 05.003.433S;  
05.003.434S; 05.003.435S; 05.003.442S; 05.003.443S; 05.003.451S; 05.003.452S; 05.003.453S;  
05.003.454S; 05.003.500S; 05.003.521S; 05.003.523S; 05.003.524S; 05.003.525S; 05.003.533S;  
05.003.534S; 05.003.535S; 05.003.553S; 05.003.554S; 05.003.703S; 05.003.721S; 05.003.723S;  
05.003.733S; 05.003.734S; 05.003.735S; 05.003.753S.

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**Appendix No. 2 – Authorized Representative**

**INSTRUCTION**

**for Use of the Medical Devices**

« Oscillating saw blade, sterile, disposable for Endoprosthetics”  
produced by Gebr. Brasseler, GmbH & Co. KG, Germany

**Appendix 2 (for Russian Federation)**

**INFORMATION ABOUT THE AUTHORIZED REPRESENTATIVE FOR MEDICAL DEVICES  
CIRCULATION IN RUSSIA:**

The limited liability company "MedTestConsulting" (LLC " MedTestConsulting")  
Address: Russia, 123181, Moscow, Kulakova str., bldg. 1, korp. 1, office 193  
Phone numbers: +7 (495) 750-29-01, +7 (903) 729-77-14  
e-mail: medtestconsult@mail.ru

This appendix is an integral part of the document "Instructions for the use of medical devices  
« Oscillating saw blade, sterile, disposable for Endoprosthetics” produced by Gebr. Brasseler  
GmbH & Co. KG, Germany

Gebr. Brasseler GmbH & Co. KG

Lemgo, 20/04/2020

Gebr. Brasseler GmbH & Co. KG  
Trophagener Weg 25 · 32657 Lemgo  
Postfach 160 · 32631 Lemgo · Germany.

*Klaus Rübesamen*

Klaus Rübesamen  
Managing Director

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  
\_\_\_\_\_  
( Nagel )

RU

;- Symbols to be used with medical  
e supplied - Part 1: General

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Manufacturer's catalogue  
identified.

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Manufacturer's batch code so that the

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Indicates a medical device that has

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5.2.8: Indicates a medical device that  
damaged or opened.

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Medical device manufacturer, as  
of 98/79/EC.

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Indicates the date when the medical device

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Indicates the date after which the medical device is

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Indicates a medical device that is intended for one  
single procedure.

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PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

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**Authorized Representative in Russia:**

LLC MedTestConsulting  
Russia, 123181, Moscow, ul Kulakova, 1-1-193  
Phone: +7 (495) 750-29-01, +7 (903) 729-77-14  
medtestconsult@mail.ru

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IFU # - Rev # ...

**(EN) Instructions for use.**

Oscillating saw blade, sterile, disposable, for  
Endoprosthetics

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EN

Symbols out of ISO 15223-1: Medical devices - Symbols for device labels, labelling and information to be used in accordance with requirements

**REF**

Catalogue number, Ref. no. 5.1.6: Indicates the number so that the medical device can be identified

**LOT**

Batch code, Ref. no. 5.1.5: Indicates the manufacturer batch or lot can be identified.

**STERILE R**

Sterilized using irradiation, Ref. no. 5.2.4: Indicates the device has been sterilized using irradiation.

Do not use if package is damaged, Ref. no. 5.2.3: Indicates the device should not be used if the package has been damaged

Manufacturer, Ref. no. 5.1.1: Indicates the manufacturer defined in EU Directives 90/385/EEC, 93/42/EEC

Date of manufacture, Ref. no. 5.1.3: Indicates the date the device was manufactured.

Use-by date, Ref. no. 5.1.4: Indicates the date the device is not to be used.

Consult instructions for use, Ref. no. 5.4.3: Indicates the user should consult the instructions for use.

Do not re-use, Ref. no. 5.4.2: Indicates the device is for single use, or for use on a single patient during a single procedure

er to consult the instructions  
s warnings and precautions  
n the medical device itself.

ity  
le by or on the order of an

lant site, work on bones or  
ure. Oscillating saw blades  
prosthetics. The  
ied medical professionals  
s conducted by duly  
ates oscillating movements.

apes of working part  
different connection types  
Ti). The thickness of the saw

paration of bones according  
intended for treatments in

olves risks outweighing the  
sed if any other  
gery in general.  
cted.

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Caution, Ref.no. 5.4.4: Indicates the need for the user to be aware of important cautionary information such as warnings that cannot, for a variety of reasons, be presented on the device.

Regulated by European law: CE marking of conformity

**Rx Only**

Non-regulated symbol:

Caution: Federal law (USA) restricts this device to sale only to a licensed physician

Non-regulated symbol: Quantity

Non-regulated symbol: Material used

**GS 1**

Barcode: GS1 format

Non-regulated symbol: Distributor

Material: Stainless steel

**Purpose:**

The instruments are used for preparation of the implant site. The instruments are used for preparation of the implant site separation of bones according to the chosen procedure. The instruments are intended for treatments in foot / knee / hip endoprosthetics. The instruments are used in operating theatres by qualified surgeons (surgeons, orthopaedists). Installation of saw blade is done by qualified medically trained users (surgical assistants).

The blades are driven by a power system that generates high speed oscillations.

The saw blades are available in different sizes and shapes (parallel or flared) with sharp cutting teeth and with different cutting angles (DePuy Synthes/Zimmer Biomet, Stryker, Aesculap 3Ti). The thickness of the blades varies between 0,89 – 1,47 mm."

**Intended use:**

Preparation of the implant site; work on bones or separation of bones according to the chosen procedure. Oscillating saw blades are intended for use in foot/knee/hip endoprosthetics.

**Contraindications:**

- Operations in which the use of a motorized saw involves a high risk of complications. A motorized saw should not be used in operations where the contraindications apply that prohibit orthopedic surgery
- The contraindications of the implant must be respected.

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**Proper use:**

- Instruments may only be used in the appropriate h with connection types: DePuy Synthes/Zimmer Biom
- Prior to use, make sure that the instrument is safe proper position.
- Contact of the saw blade with metal and other met avoided.
- The saw blade should be actuated in / at the saw g made with the bone.
- The thickness of the saw blade must be inferior to 1 template.
- Free mobility of the saw blade in/at the template n
- The instrument should be carefully inspected before functionality.
- Single-use only. Do not re-use. Instruments must be operation.
- The instructions for use for the power system used
- In case of breakage of the metallic instruments frag can be located by medical imaging procedures.

**Packaging:**

- The present instrument is sterile packed.
- Sterilization was placed by means of gamma irradiat
- The instruments need to be discarded if the expiry c packaging is opened or damaged or in case of illegibl
- The patient labels must be added to the patient's he etc. after initial opening.

**Storage, transportation and usage:**

Prior to use, instruments must be stored in their origi from direct sunlight, dust, moisture, sources of reco temperature .

**Operating conditions:**

Temperature: +5°C to +45°C

Relative humidity: 30% to 80%

The instruments should be stored in a dry state at a t +50°C and humidity of 10 to 95%.

**Transportation conditions:**

The instruments are transported by all modes of trans accordance with the rules of transportation applicable at a temperature from -50°C to +50°C and humidity fr

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**Warnings:**

- Avoid jamming or bending of the instrument. Do not use the instrument if jamming or bending leads to excessive heating, jamming of the saw blade.
- Do not allow the saw blade teeth to come in contact with the skin or other metal objects while performing oscillating motion. This leads to damage to the sawblade or tools as well as to the patient, especially to bones and surrounding tissue (incl. nerves and blood vessels). This increases the risk of inadvertent creation of metal chips.
- Excessive contact pressure must be avoided. This leads to overheating of tooth surfaces, reduced service life of instruments, necrosis, fracture of the blade, jamming of the blade, and an undesirable rough and uneven surface.
- Due to wear of the material, a safe and risk-free re-use is not guaranteed. Reprocessing of single-use device may lead to contamination and can result in infection and irreversible damage.
- Insufficient cooling leads to clogging of the teeth and necrosis, reduction of the service life.
- Shifting of the template or slipping of the saw blade may lead to mistakes.
- Damaged instruments (scratches, dents, blunt blades) may result in the fracture of the instrument or treatment of the patient getting injured by metal fragments left behind in the tissue.

**Precaution:**

- Single-use only. Do not re-use.
- Do not use the instrument if damaged or in case of corrosion or illegible markings.
- Do not use the instrument if the expiry date indicated on the packaging has been exceeded, the packaging is opened or damaged, or if there are illegible markings or corrosion.
- Sufficient cooling (saline solution) should be provided during use.

The instruments are not provided with a measuring device, an internal energy source and do not incorporate soft tissue. They do not contain any medicinal product constituents, ingredients or derivatives made of human blood. They do not emit (ionizing) radiation or create magnetic fields. As a result, no special precautions are expected.

**Environmental risks of medical devices:**

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**List of standards:**

DIN EN ISO 13485, DIN EN ISO 10993, DIN EN ISO 13485, DIN EN ISO 14971

**Notification of unwanted events/occurrences:**

All serious incidents occurring in connection with the use of the instrument should be notified to the manufacturer and to the competent national authority.

occurs in accordance with  
nPiN 2.1.7.2790-10).  
pointed and sharp objects

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to use to ensure that it is  
ible for the use of the  
: user, Gebr. Brasseler  
images, particularly if these  
or use or warnings as well

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**Disposal:**

Disposal of instruments that have become unusable or the domestic legislation of the Russian Federation (San Medical waste hazard class B. Potentially infectious, po must be disposed off in disposable containers.

**Quality related issues** shall be addressed to:

LLC MedTestConsulting  
Russia, 123181, Moscow, ul Kulakova, 1-1-193  
Phone: +7 (495) 750-29-01, +7 (903) 729-77-14  
medtestconsult@mail.ru

**Liability:**

The user is responsible for checking the product prior suitable for the intended purpose. The user is responsible for instruments. In case of contributory negligence by the partially or totally declines liability for all resulting damages are due to non-observance of our recommendations if as inadvertent misuse by the user.

**Trademark Reference:**

The third party trademarks used herein are trademarks of their respective owners.

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© Gebr. Brasseler GmbH & Co. KG