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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LIFETECH SCIENTIFIC (SHENZHEN) CO., LTD. on the annexed DECLARATION is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Hu Na

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(Date: Nov. 13, 2018)

Statement

November 9, 2018

We

Lifetech Scientific (Shenzhen) Co., Ltd

Floor 1-5, Cybio Electronic Building, Langshan 2nd Street,

North Area of High-tech Park, Nanshan District,

Shenzhen, 518057

P.R. China

Here state that:

The stamp in the following page is legally valid and only used for product registration in Russia and not for other use.

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Lifetech Scientific (Shenzhen) Co., Ltd.,
Floor 1-5, Cybio Electronic Building,
Langshan 2nd Street, North Area
of High-tech Park, Nanshan District,
Shenzhen, 518057, P.R. China

LiLing/ Regulatory
Affairs Supervisor

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CeraFlex¹ ASD Closure System, in different versions

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1 Device Description

The CeraFlex™ Atrial Septal Defect (ASD) Closure Systems include CeraFlex™ ASD Occluders and the SteerEase™ Introducer. The CeraFlex™ ASD Occluders have two types of occluders, CeraFlex™ ASD Occluder and CeraFlex™ Multi-fenestrated (Cribiform) ASD Occluder as shown in **Figure 1**. The CeraFlex™ ASD Occluders are percutaneous, transcatheter closure devices for the non-surgical closure of atrial septal defects. The CeraFlex™ ASD Occluder is self-expandable double disc device, which is made of a nitinol wire mesh that is shaped into two flat discs and a waist between the two discs. PET membranes sewn into each disc and the waist (no PET membrane sewn into the waist for multi-fenestrated ASD occluder) help to seal the hole and provide a foundation for growth of tissue over the occluder after placement. All metallic structures are plated with titanium nitride (TiN) to improve the device's biocompatibility.



Figure 1. The Structure of the CeraFlex™ ASD Occluders

The CeraFlex™ ASD Occluders must be used in combination with the SteerEase™ Introducer. The introducer contains a coil reinforced sheath, a dilator, a loader, a haemostatic valve and a delivery cable. It is used to advance the CeraFlex™ ASD Occluder to the proper position. When the occluder is released from the sheath, the disc expands on each side of the defect. At the same time the waist expands to its pre-determined diameter and closes the hole between the left and right atria. The CeraFlex™ ASD Closure Systems consist of an ASD occluder, a Delivery cable, a Loader, a Haemostatic valve, a Sheath and a Dilator. The system (Refer to **Figure 1** and **Figure 3**) is comprised of

a) Left atrial disc view of CeraFlex™ ASD Occluder

b) Right atrial disc view of CeraFlex™ ASD Occluder

Figure 2. The Structure of the CeraFlex™ ASD Occluders

- **ASD Occluder:** The device for closing ASD. The CeraFlex™ ASD Occluder is connected with a delivery cable by a loop connection through holes in the proximal portion of the occluder, adjacent to the proximal hub. The CeraFlex™ ASD Occluder is preassembled with delivery cable.
- **Sheath:** The sheath is used to advance the device to its desired position within the heart.
- **Dilator:** The dilator is used to ease penetration of tissues and the vessel wall.
- **Hemostatic valve:** The haemostatic valve at the proximal end of the sheath minimizes bleeding. The sideport with the flexible extension tube and stopcock is used to flush the system.
- **Loader:** The loader is used to introduce the occluder with the attached delivery cable into the sheath.

Delivery cable with Plastic vise/Handle (proximal): The delivery cable is used to advance (push) the occluder through the sheath, holding it in position while the sheath is pulled back to deploy the occluder. The Handle and delivery cable are also used to retrieve and/or reposition the occluder if the size, position or the expansion of the occluder is deemed to be unsatisfactory. The handle is used to facilitate direction control and serve as "release control handle" for disconnecting (releasing) the occluder from the delivery cable.


Figure 3. CeraFlex™ ASD Occluder System and the SteerEase™ Introducer

2 Indications and Usage

The CeraFlex™ ASD Occluder is a percutaneous, transcatheter, closure device intended for the occlusion of atrial septal defects (ASDs) and multi-fenestrated atrial septal defects in secundum position or patients who have undergone a fenestrated Fontan procedure and who now require closure of the fenestration. The CeraFlex™ ASD Occluder applies to closing secundum atrial septal defects, and the CeraFlex™ Multi-Fenestrated (Cribriform) ASD Occluder applies to closing multi-fenestrated atrial septal defects. Patients indicated for ASD closure have echocardiographic evidence of (fenestrated) ostium secundum atrial septal defect and clinical evidence of right ventricular volume overload.

3 Contraindications

- Patients who have extensive congenital cardiac anomaly which can only be adequately repaired by cardiac surgery.
- Patients who have sepsis within one month prior to implantation, or any systemic infection that cannot be successfully treated prior to device placement.
- Patients known to have a bleeding disorder, untreated ulcer or any other contraindications to aspirin therapy unless another anti-platelet agent can be administered for 6 months.
- Patients who have demonstrated intra-cardiac thrombi on echocardiography (especially left atrial or left atrial appendage thrombi).
- Patients whose size (i.e., too small to tolerate TEE probe, catheter size, etc) or condition (active infection, etc.) would cause the patient to be a poor candidate for cardiac catheterization.
- Patients where the margins of the defect are less than 5mm from the coronary sinus, AV valves or right upper lobe pulmonary vein (for CeraFlex™ ASD Occluder).
- Any patients for whom the radius of the device is greater than the distance from the central defect to the aortic root or superior vena cava (for the CeraFlex™ Multi-Fenestrated (Cribriform) ASD Occluder)

4 Warnings

- Patients allergic to nickel may suffer an allergic reaction to the device.
- The CeraFlex™ ASD Closure Systems must only be used in hospitals, which have experience with invasive treatment of congenital heart disease, and used by physicians who are trained in transcatheter defect closure techniques.
- Physicians must be prepared to deal with emergency situations which require removal of embolized devices that result in critical hemodynamic compromise. This includes the availability of an onsite surgeon.
- Embolized devices must be removed. Embolized devices should not be withdrawn through intracardiac structures unless they have been adequately collapsed within a sheath.
- The CeraFlex™ ASD Closure Systems are for single use only. Do not reuse or resterilize. Structural integrity and/or function may be impaired or lost through cleaning, resterilization, or reuse and may cause adverse reactions. Lifetech will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the ASD closure systems.
- Do not release the CeraFlex™ ASD Occluder from the delivery cable if the device does not conform to its original configuration or if the device position is unstable. Recapture the device and redeploy. If still unsatisfactory, recapture the device and replace with a new device.
- Implantation of this device may not eliminate the need for Coumadin in patients with ASD and paradoxical embolism.
- Implantation requires the use of echocardiographic imaging (TTE, TEE or ICE). Do not inflate the sizing balloon, if used, beyond the cessation of the shunt (i.e., stop-flow) or the visualization of a small waist in the balloon. Do not overinflate. (Balloon Sizing may not be suitable for the CeraFlex™ Multi-fenestrated (Cribriform) ASD Occluder.)
- Do not select a device which size is 1.5 times larger than the ASD diameter measured by the echocardiographic – device prior to balloon sizing (for CeraFlex™ ASD Occluder).
- The Occluder is made of NiTi wires, which is a material sensitive to temperature. Therefore, some occluders might present malformation at low temperature. Keeping the temperature of saline water and operating room above 24 °C would help occluders to exhibit the desired shape.

A. Cause for peanut malformation:

- a) Over-sized Occluder is selected.
 - b) Retract and release the occluder from loader or sheath repeatedly.
 - c) The temperature of saline water or the operating room is too low when an occluder is loaded into the loader. This may result in peanut malformation of the occluder, which cannot recover immediately after the occluder is implanted.
- B. Cause for tulip malformation:**
- a) When a force is applied on the left disc or when the left disc is being released against the left atrial wall, the top of the left disc might be distorted resulting in tulip malformation.
 - b) During retraction of the occluder, if the left disc of the occluder is pushed against the left atrial wall or the wall of saline water container, tulip malformation may occur when the left disc is released afterwards.
- C. How to prevent malformation:**
- a) Some occluders might present malformation at low temperature. Keeping the temperature in the operating room and the saline water above 24°C would help occluders to exhibit the desired shape.
 - b) Pull the left disc of the occluder axially to extend the occluder, and then retrieve the occluder into the loader.
 - c) Do not retrieve the occluder into loader against the wall of saline water container.
 - d) Do not retract and release the occluder from loader or sheath repeatedly.
 - e) Do not release the occluder against the left atrial wall.
 - f) The membrane might be jammed when retrieving the occluder into the loader, which is another possible reason for malformation. Pay attention to the condition of the membrane when retracting the occluder into loader.
- D. How to correct the device, when malformation is occurs**
- a) If malformation occurs at the left disc, retrieve the occluder, then push the occluder completely out of the sheath in the left atrium, use the tip of sheath to push occluder until the left disc is rested against the left atrial wall to correct the shape, and then release the occluder in the defect area.
 - b) If malformation occurs at the right disc, advance the sheath and gently push the tip of sheath against the right disc to correct the shape.
 - c) If a peanut or tulip malformation is observed after the occluder is released, retrieve the occluder and immerse it in warm water (above 24°C). Extend and retract the occluder repeatedly to relax the NiTi wires until the shape of the occluder is recovers.

5 Precautions

- The use of this device in patients with patent foramen ovale has not been studied.
- The CeraFlex™ ASD Closure Systems are for single use only. Do not reuse or resterilize.
- Patient Selection

Certain patients may be at higher risk for complications such as tissue erosion and device embolization. If higher risk patients have devices implanted, closer follow-up is warranted. Higher risk patients include the following:

 - Patients with deformation of the device at the aortic root
 - Patients with high defects (minimal aortic and superior rims)
 - Patients with IVC rim deficiency (risk of device embolization)
- Pre-procedure, procedure and post-procedure
 - Pre-procedure

Aspirin (3-5 mg/kg/day) is to be started at least 24 hours prior to the procedure. In the rare case of aspirin intolerance, two times 200 mg of Ticlopidin are given. Cephalosporin therapy is optional.
 - Procedure

Patient should be fully heparinized throughout the procedure with a minimum active clotting time (ACT) of 250 seconds prior to device insertion.

Transesophageal echocardiography (TEE) or similar imaging equipment is recommended as an aid in placing the ASD Occluders in the proper position. If TEE used, the patient's esophageal anatomy must be adequate for placement and manipulation of the TEE probe
 - Post-procedure

Patients are required to take appropriate endocarditis prophylaxis for 6 months following device implantation. The decision to continue endocarditis prophylaxis beyond 6 months is at the discretion of the physician.

Patients, who had any residual shunt, should be inspected by echocardiographic examination for the residual shunt until the defect closed completely.
- MR conditional

A patient with an implanted ASD Occluders can be scanned safely immediately after placement of the device under the following conditions:

 - Static magnetic field of 3 T or less
 - Spatial gradient magnetic field of 720 G/cm or less
 - Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning

Note: MRI image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

6 Potential Adverse Events

Placement of the CeraFlex™ ASD Occluder involves using standard interventional cardiac catheterization techniques. The following adverse events might occur from interventional cardiac catheterization techniques or from device.

- ◇ Allergic reaction
- ◇ Arterial fistula or left or right atrial fistula
- ◇ Septal tearing
- ◇ Injury to the nerve or vessel wall
- ◇ Incomplete sealing of the defect
- ◇ Mitral insufficiency
- ◇ Transient ischemic attack
- ◇ Fever
- ◇ Arrhythmias
- ◇ Thrombus formation
- ◇ Device migration
- ◇ Cardiac perforation
- ◇ Infection
- ◇ Cardiac tamponade
- ◇ Valvular regurgitation

7 Product Features

The CeraFlex™ ASD Occluders and the introducers are available with specifications in the following (Table 1).

Figure 4. CeraFlex™ ASD Occluders

Table 1. List of specifications of CeraFlex™ the ASD Closure System

Occluders specification	ØD2	ØD1	ØD	L	Recommended introducers
CeraFlex™ ASD Occluder					
LT-ASDf-06	18	14	6	4	SFA8F-f
LT-ASDf-08	20	16	8	4	SFA8F-f
LT-ASDf-10	22	18	10	4	SFA8F-f
LT-ASDf-12	26	22	12	4	SFA9F-f
LT-ASDf-14	28	24	14	4	SFA10F-f
LT-ASDf-16	30	26	16	4	SFA10F-f
LT-ASDf-18	32	28	18	4	SFA10F-f
LT-ASDf-20	34	30	20	4	SFA12F-f
LT-ASDf-22	36	32	22	4	SFA12F-f
LT-ASDf-24	38	34	24	4	SFA12F-f
LT-ASDf-26	40	36	26	4	SFA12F-f
LT-ASDf-28	42	38	28	4	SFA12F-f
LT-ASDf-30	44	40	30	4	SFA14F-f
LT-ASDf-32	46	42	32	4	SFA14F-f
LT-ASDf-34	50	44	34	4	SFA14F-f
LT-ASDf-36	52	46	36	4	SFA14F-f
LT-ASDf-38	54	48	38	4	SFA14F-f
LT-ASDf-40	56	50	40	4	SFA14F-f
LT-ASDf-42	58	52	42	4	SFA14F-f
CeraFlex™ Multi-Fenestrated ASD Occluder					
LT-ASDf-MF-1818	18	18		3	SFA9F-f
LT-ASDf-MF-2518	18	25		3	SFA10F-f
LT-ASDf-MF-2525	25	25		3	SFA10F-f
LT-ASDf-MF-3030	30	30		3	SFA12F-f
LT-ASDf-MF-3525	25	35		3	SFA14F-f
LT-ASDf-MF-3535	35	35		3	SFA14F-f
LT-ASDf-MF-4040	40	40		3	SFA14F-f

8 Individualization of Treatment

● Patients Selection

Device placement should only be attempted in those patients with sufficient rim around the defect to allow stable seating of the device.

● Patients with multiple ASDs

- Closure of multiple ASDs should only be attempted by physicians who have sufficient experience (more than 10-15 cases) to undertake more technically challenging procedures.

4 English

- When two large ASDs separated by more than 7 mm rim of tissue, then implantation of two devices may be justified.
- When there are multiple ASDs close to each other, one device may be used to cover all defects by placing the device in the largest defect
- Device placement and size selection
 - Device placement is required to be done with the assistance of TEE or similar imaging equipment.
 - Select an occluder about 2-4 mm in diameter larger than the ASD defect (refer to Table 1).
- Use in specific populations
 - Pregnancy – Care must be taken to minimize the radiation exposure to the fetus and the mother.
 - Nursing mothers – Although appropriate biocompatibility testing has been conducted for this implant device, there has been no quantitative assessment of the presence of leachables in breast milk.

9 Directions for Use

- Perform general or local anesthesia. If TEE used, the patient has to be intubated endotracheally.
- Following percutaneous puncture of the femoral vein, perform a standard right heart catheterization.
- Perform an angiogram in order to demonstrate the atrial communication. Catheterize the left atrium using a 45° LAO position and a 35-45° cranial angulation; inject contrast medium into the left atrium or into the right superior pulmonary vein.
- Introduce a 0.038-inch exchange "J" tip guidewire into the left atrium. If desired, insert a compliant sizing balloon catheter over the exchange wire into the left atrium and measure the defect size.
- Sizing the defect: measure the diameter of the ASD by using a compliant sizing balloon, which is filled with diluted contrast, placed across the ASD, and inflated until a waist in the balloon is visible on fluoroscopic imaging. The balloon is gradually inflated until shunting is eliminated on the Doppler evaluation.
- For multi-fenestrated ASD, device size and placement (the fenestration in which the device will be placed) is based on the location of the fenestration. The use of transesophageal echocardiography or intracardiac echocardiography is mandatory in positioning the sheath in the most centrally-located fenestration. Once the most centrally-located defect has been crossed, measurements should be made to determine the distance to the outer rim of the farthest fenestration, and the device size selected should be 2:1 of this measurement. The device selected must completely cover all fenestrations.
- Once the diameter of the defect has been determined, select an occlusion device equal to or, if the confirmed size is not available, one size larger than the defect. Select the proper size sheath to accommodate the occluder.
- Remove the balloon catheter leaving the exchange guidewire in place.
- Flush the dilator and sheath with saline solution, and then insert the dilator into the sheath, tighten the clamp nut to ensure the dilator fasten with sheath securely. Advance the sheath with dilator over the exchange guidewire through the septal defect into the left superior pulmonary vein. Remove the exchange guidewire and dilator (Refer to Figure 5).
- Clockwise screw the loader into the haemostatic valve until it stops, and then check connection is reliable.
- Infuse the saline solution from the stopcock to flush the haemostatic valve and the loader, while keeping the tip of the loader immersed under the saline solution.
- Carefully check the connection between the occluder and the delivery cable (Refer to Figure 6). Immerse the device and loader in saline solution and pull the device into the loader with a quick jerking motion (Refer to Figure 7)
- Twist the seal cover of the haemostatic valve to suitable degree (Refer to Figure 8), and then repeatedly push and pull to make sure the tension is appropriate.

Warning:

If the seal cover is too tight, then the cable will be hard to push. However, if the seal cover is too loose, there will be a risk of aeroembolism.

- Flush the device via the side arm.

Figure 5. Illustration of delivering sheath into heart


Figure 7. Illustration of loading the device into a loader


Figure 6. Illustration of CeraFlex™ ASD Occluder System

Figure 8. Illustration of securing the seal cover of the haemostatic valve

- Keep the sterile heparinized saline outflowing from the loader to eliminate air in the loader. Simultaneously, insert the loader into the haemostatic valve of the sheath gently until it stops, and then tighten the clamp nut of loader. Advance the occluder into the sheath by pushing (not rotating) the delivery cable. **Do not remove the loader from the sheath when advancing the device.**

Warning:

If there is air in the loader and sheath when advancing the device, it will result in an air embolism. If the loader is removed from the sheath when advancing the device, there will also be risk of an air embolism.

- Under fluoroscopic and echocardiographic guidance, deploy the LA disc and waist and pull the device gently against the atrial septum, which maybe felt and also observed by ultrasound imaging. Then immobilize the delivery cable, pull the sheath back and deploy the right atrial disk. Pull the sheath back by approximately 15 cm. A gentle "forward and backward" motion with the delivery cable verifies a secure position across the atrial septal defect, which can also be observed by TTE or TEE (Refer to Figure 10).
- Confirm the device position and evaluate residual shunts. Perform angiography or echocardiography. If the position is unsatisfactory or there is obvious residual shunt, stabilize the delivery cable and re-advance the sheath until the device is completely within the sheath. Reposition the device and deploy it again, or remove the device from the patient. If device malformation has occurred, deploy two discs in the left atrium and it may be necessary to push gently against the atrial wall to recover the original device configuration.

Warning:

When withdrawing the occluder back into the sheath, pull back on the delivery cable, and do not push the sheath. Otherwise, the tip of the occluder maybe stuck outside the sheath, and then the occluder cannot be resheathed.

- Release the device when the position is satisfactory. To release the device, press the blue button and move back, while maintaining the orange button press state.
- Hold the sheath to touch the hub of the occluder, and then pull back the cable gently until it separates from the occluder. Finally remove the delivery cable and sheath from the patient (Refer to Figure 9).
- Use TTE or TEE to confirm device placement, assess for residual shunts, obstruction, or regurgitation induced by the device. Perform an angiogram to assess for residual flow through the device.


Figure 9. Illustration of releasing the occluder

1. Press the lock button as the arrow of ① (the color of the button is orange).
2. Move the control button as the arrow of ② (the color of the button is blue).

Figure 10. Illustration of implanting the occluder

10 Post-procedure

- All patients should be kept overnight for observation.
- Anticoagulation treatment with heparin for 24 hours.
- Oral anticoagulation with Aspirin 3-5 mg/(kg.d) for children, and 3 mg/(kg.d) for adults, 6 months; conditionally, plus Plavix 75 mg/d for some adult patients.
- Antibiotic therapy
- Re-examination by transthoracic echocardiography (TTE) at 24 hours prior to discharge.
- Patients with high risk should be followed more closely, including the following:
 - Clinical follow-up with echocardiogram one week following device implantation.
 - Education of patients about the higher risk and the need for echocardiography in case of symptoms (i.e., chest pain or shortness of breath)

11 Expiry Date

The ASD Closure Systems are sterilized by ethylene oxide. The shelf life is indicated on the labels, and the manufacturing date and use by date are marked on the label. **Do not** use an expired product.

12 Package and Label

The CeraFlex™ ASD Closure Systems are supplied sterile. The CeraFlex™ ASD Occluder is pre-connected with delivery cable by loop connection and pre-assembled with its loader and Haemostatic valve. The sheath and dilator are packaged separately. The package content is shown in Table 2.

Table 2. The packaging contents of the ASD Closure Systems

Package	CeraFlex™ ASD Closure Systems	
	CeraFlex™ ASD Occluders	SteerEase™ Introducers
Contents	an ASD Occluder, a Delivery Cable, a Loader, and a Haemostatic Valve	a Sheath and a Dilator

The ASD device is fixed and protected on a PETG tray, then sealed into inner and outer Tyvek/filmpouches, where a label and a sterilization indicator are attached. The product is sterilized by ethylene oxide gas and put in a box with the Instructions for Use, a patient card, a customer feedback form, and a certificate of compliance. An outer label is attached to the box.

The sheath and dilator of the SteerEase™ Introducer are fixed and protected on a scale board, then sealed in two dialyzing pouches, where a label and a sterilization indicator are attached. The product is sterilized and put in a box with the Instructions for Use, a customer feedback form, and a certificate of compliance. A label is attached to the box.

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LIFETECH SCIENTIFIC (SHENZHEN) CO., LTD. on the annexed STATEMENT is genuine.

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17 January, 2019

We

Lifetech Scientific (Shenzhen) Co., Ltd
Floor 1-5, Cybio Electronic Building, Langshan 2nd Street,
North Area of High-tech Park, Nanshan District,
Shenzhen, 518057
P.R. China

Here state that:


The document only used for product registration in Russia and not for other use. All the contents are true and lawful.

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лица/name of responsible person)

Ling Li/Regulatory Affairs Supervisor

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(Lifetech Scientific India Private Limited)